VOCO

Meron
GB Instructions for use **CE**⁰⁴⁸²
in accordance with EN ISO 9917

Meron is a glass ionomer cement (glass polyalkenoate cement) for luting.

Meron has a low solubility in the mouth and a low acidity.

Fields of application:

For attaching crowns, bridges, inlays, onlays, pins, posts and orthodontic bands.

Method of use:

Clean enamel, dentine and metal surfaces thoroughly and dry completely. Make sure not to leave any traces of other materials which could impair the properties of **Meron**. Temporary work should be attached with an eugenol-free cement. Before attachment clean the core, rinse (e.g. with H₂O₂) and dry.

When using **Meron** for the first time, remove cap from bottle and fit the pipette attachment supplied.


Shake the powder well before use!

The working temperature should be 15°C - 23°C.

Use one drop of liquid for one level measuring scoop of powder (powder/liquid ratio = 3,0 : 1 g/g). Put powder and liquid on a suitable glass plate or a suitable mixing pad.

When dispensing the liquid, hold the dropper bottle vertically.

It is recommended to divide the powder in 3 portions and to mix it portion by portion into the liquid, using a solid plastic spatula or another abrasion-resistant spatula.

Application:

Meron has a working time of approx. 3 min and sets after 5 - 7 min (3 - 5 min in the mouth). Higher temperatures shorten the working time, low temperatures prolong it.

The correct consistency is achieved when **Meron** is applicable to the core with a hard brush. It is essential to avoid contact with water and saliva during application and setting.

Indications and Precautions:

Excess material can be removed after 10 - 15 min, preferably with an explorer. Desensitize strongly ground and sensitive cores prior to a temporary attachment. Remaining varnish and other substances (e.g. isolating agent) can affect the adhesion and should be removed completely before using **Meron**.

Do not apply the cement excessivley and maintain sufficient dentine thickness to avoid pulp irritations. If necessary, pulp protection should be applied in the proximity of the pulp or a core build-up should be carried out, if indicated.

Rinse instruments with cold water immediately after use and before the setting of the cement residues. After removal of orthodontic bands, cement residues can be removed gently and rapidly with an ultrasonic scaler.

Meron contains polyacrylic acid, fluoro-silicate and parabens. Do not use **Meron** in cases of hypersensitivities against these ingredients.

Storage information:

Close powder container immediately after use, in order to avoid infiltration of humidity which may result in a loss of compressive strength. Store at temperatures 4°C - 23°C.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VOCO GmbH Phone +49 (4721) 719-0