# **DENTA**FLUX

## **THE MACHINE**


(Adaptaciones termoplásticas)

MANUAL DE INSTRUCCIONES

INSTRUCTION MANUAL

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[www.dentaflux.com](http://www.dentaflux.com)

J Ripoll S.L.  
Río Tietar 20 Pol Nogal  
28110 Algete Madrid

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
Please read carefully this instructions before connecting to the Light source.

The use of the machine is easy and can be learned in few minutes with this instructions.

- 1.-heater
- 2.-heater andel
- 3.-axis vertica
- 4.-upper frame
- 5.-lower frame
- 6.-frame closer
- 7.-suction base
- 8.-frame andel
- 9.-body (motor)
- 10- heater switch
- 11- heater lamp
- 12-suction motor switch

## INDEX

Important notes

Instalation

Basic instructions

Uses:

Baseplate

Custom trays

Other rigid uses

Surgery trays

Bruxism flexible trays and sport protectors

Temporary crowns and bridges

Adaptations for surgery procedures

Duplicates

Whitening trays

Sheets avaiable

## IMPORTANT NOTES

- 1.- When using mouthguard materials it is recommended to use a silicon spray over the model and over the frame's rubber. It make longer self life of the rubber and make easier the separation between the plastic and the model. Anyway the experience Hill help you.
- 2.- With the time the rubber of the frame should be changed. Remove also the adhesive rest. Any abrasive can help you. To join the new rubber use an adhesive apropiate. You can get it in the depot you purchased the machina.
- 3.- The clear materials (rigid and transparent) can absorbe ambiantal humidity, This create bubbles when heated. If it is not too wetted, can be dried in a home use oven keeping at 60-70°C during 1 hour, as it is said in the plastic sheets box.
- 4.- Dentaflux the Machine just use high quality materials, but anyway after years of use the heater or suction motor could need a change. The spare pieces are aviable for an inmediate repair.
- 5.- The vertical axis should be oiled frequently with an Light oil.
- 6.- If you have any doubt that your depot cannot solve you can freely contact us at 34916289399 (english, french or spanish)

## INSTALATION

A

The machine is ready to use . Just few steps are necesry.


B

Get the handles from the bag, **place two in the back zone of the frame** (near the vertical axis) **and one in the front part**, in the right or left hole depending you are lefty or not. There are another handel to be placed in the heater head.

C


Connect to the Light source

D


Ones the heater works just touch the head with the handle. This piece gets very hot.

## BASIC INSTRUCTIONS


1.- Turn the head to the back and switch on the heater. The unit will need a couple of minutes to heat. Meanwhile you can do the steps 2 to 5. The lamp indicates if the heater is connected or not.

2.- Place the selected plastic sheet inside the frame and close it. To open the frame you should free it from frame closer.



3.- **Place two handles near the vertical axis and just one to open and close the frame.** By the two handles near the vertical axis, put the frame up until a click is heard. You will see how it gets fixed. Sometimes the

new machines can be hard to elevate the frame. Please check the screws, they could be too tightened.

4.- Place the model (perhaps two models could fit) over the suction base. The machine let two models if they are small and should be placed in diagonal with a small separation between them. The models must be always trimmed leaving just the area to duplicate. When using rigid sheets, the retentions must be covered with refractory mass to let separate the model from the plastic after the adaptation is done. If a hole is done in the palatal area, it helps the suction to a better work.


5.- Place the heater over the plastic.

**PRECAUTION:** the heater can be turned only with the handle to avoid harm.

Do not let the plastic heating unattended, if overheated the plastic will be useless and with time it will melt, goes through the holes towards the motor and destroy it.

6.- When the plastic starts to fall, switch on the suction, and by the two posterior handles set down the frame. The

experience will show when the plastic is ready for adaptation. The thinnest ones (0.20 or 0.5mm) just need to be soft if touched below. As the plastic is thicker it needs to fall. The thickest as the mouthguard 150 (3.5mm) can drop 2-3cm before adapting.

7.- LEAVE THE HEATER AND SUCTION CONNECTED DURING 10-15SG, UNTIL THE ADAPTATION IS COMPLETED.

8.- Turn back the head. Leave the plastic to cool with the suction connected during 30 sg.

9.- Switch off the motor. Take out the adaptation and leave the frame open to get cool.

10.- To separate the plastic from the model. If a mouthguard material has been used, just let it cool and separate, the plastic will keep the form, cut with scissors and finish. If a thin plastic has been used a knife can be useful to cut it. If a rigid plastic has been used, a refractory paste should be used to let separate the model from the plastic.

11.- There is a special disc to cut plastics. Depending of the sheet, a knife or a bur can be useful.

## BASEPLATES

Plastic: 1.5mm base material (.060 white), 2.0mm base material (080 white), 2.5mm base material (100 pink)

The basic technique is the common one and has been explained. Here there are just additional notes for this specific use.

### Material

There are different thickness depending of the rigidity desired in the adaptation. For lower should be used the thickers. When the plastic is removed from the model, you can see the cutting and finishing is easy, and it is not as fragil as the old baseplates. When this material is cutted with a bur, it melts around the bur. There are specials discs designhed for this use. If a specific area must be readapted, it can be done heating with a flame and adapting with the fingers.. Stabilizers are not needed. When the adaptation is cool it is stable. The sheet does not bond the model neither melt during the adaptation process. The adaptation fits so well to the soft tissue and powders are not needed.

### Refractory mass:

For all models with retentive areas, this should be covered with refractory mass. If the mass is not used will be not possible to separate the adapted plastic to the model. Wax should not be used because it will melt and go to the suction motor.

### Reinforze

Sometimes it is necessary to place a reinforze. Use the regular wire that would be fixed to the plastic in 2 or 3 points with cold cure acrylic or wax. Apply the monomer with a brush over the plastic in the lingual area and add the powder to get an additional strength.

### Finishing

The borders can be finished with pomez and a brush. In special cases the surface can be sand blastered to get more rugosity.

### To maintain a chamber:

Place a chamber over the model and bond it before adapt the plastic. With this we avoid removing the palatal area when the product comes back for finishing. We will get a uniform thickness in the palatal area. For finishing, just place a sheet wax for rugosity and continue with the common way.

## CUSTOM TRAYS


Plastic: 100 pink (2.5mm),125 blue (3.0mm) 150 blue (4mm)

THE MACHINE lets you get quickly and accurated custom trays with a low cost. It is not necessary any painting, hand adaptations, wait for polimerization or polish. The tray is just the adaptation cut.

1.- The model must be trimmed in the usual way. We recommend to trim the model as soon as possible leaving just the area to duplicate. This will easier the work.

2.- With a rectangle of any old sheet we will make the handle: fix the rectangle **vertically**, in the oclusal area, with wax. If the model is wet, can be necessary to make a little hole where is some inserted the handle, and fix it, avoiding that the piece get down when the sheet is adapted. In quick trays, the handle can be bonded to the lingual margin of the anterior teeth.

The handle can be made also with autocure resin.

3.- To leave space between the the adaptation and the model, wet foam can be used. This creates also rugosity. And in any special area to cover use refractory mass. Break the foam and pass the piece that will make the handle trough the foam.

4.- In deep angles it is recommended to cover with refractory mass if not necessary a good adaptation.

5.- Adapt the sheet. The “handle” will be fushed with the sheet, vertically. As it is hot, bend it as desired.

6.- To cut the material a disc can be used or a bur, but with with wide ways to avoid the plastic block it. DENTAFLUX has a special disc for this use. A soldering pencil is also suitable (available in any hardware shop), that will melt the plastic for cutting.

7.- Any finishing bur leaves the borders smooth. NOTE: if the tray will be used with alginate, retentions can be made with a round n°8 bur. Inside could be sanded and can be placed cords with wax

### **TECHNIQUE FOR THE FOLLOWING USES**

#### **BITE ELEVATOR** **RIGID ADAPTATION FOR BRUXIST** **SPACE MAINTAINERS** **ADAPTATION FOR PERIODONTAL STABILIZATION** **OCCLUSION CORRECTORS** **TEMPORARY PARTIAL**

Plastics: the sheets used are the clear ones. They are available in the following thickness 1.5mm (060), 2.0mm (080), 2.5mm (100), 3.0mm (120) and 4.0mm (150)

1.- The basic technique is the same but in this uses the model should be trimmed or covered with refractory mass leaving free only the area to adapt. Covering the other areas, the separation from the model will be easier.

2.- After adapt the plastic, cut the plastic until the desired level. A small sulcus can be made on the model. This will create a mark on the plastic to follow. For bite elevators, or to modificate oclusionation planes, you can leave the posterior zone longer to increase the stability. If you need to elevate the bite more than 1mm, that is the regular thickness obtained, autocure resin can be placed. Also a second adaptation can be made, after finish the first, separate from the model and make another adaptation. Finish the borders of both. Place one over the other and bond with autocure resin with the height desired and with a natural occlusion.

In temporary partials, if there are missing teeth, sintetic teeth can be placed on the model with wax, plaster or resin. For ortho retainers and space maintainers the model must be covered in all sharp areas.

NOTE: The 060, 080 or 100 pink can be used for extra hardness. The adaptation can be inserted to cure auto resin at 50°C.

### **SURGERY ADAPTATIONS**

The machine is so flexible and can make manu others uses in modern odontology. With the plastics 060 Clear (1.5mm) or 080 Clear (2.0mm) transparent adaptations can be done.

Obtain the first sample and boil it. Get an impression of the model with alginate.

Get the plaster model and adapt the sheet. Finish the adaptation in the usual way. Pomez can be useful.

NOTE Over clear materials: can absorb humidity from air. This produces small bubbles when heated. If it is not too wetted can be dried placing at 60°C 1 hour, and leave to cool.

## **SOFT BRUXISM ADAPTATIONS. SPORT PROTECTORS**

Plastic: 120 mouthguard (3.0mm), 150 mouthguard (3.5mm). Also coloured.

1.- Trim the model to leave it as low as possible.

NOTE: If you want to mark the adaptation, for example into a sport team, marks can be done on the model in the lingual area near the gingival margin. This will be duplicated in the protector.

2.- You can use silicone spray over the model. This produce more transparent adaptations.

3.- Adapt the sheet

4.- Can be cooled quickly under cool water.

5.- Cut the material with hard scissors.

6.- The borders can be finished with any usual piece.

7.- If you want to register occlusion areas, the adaptation should be mounted on the articulator, heat the area to soften it, and close the articulator marking the plastic until the height desired.

## **TEMPORARY CROWNS AND BRIDGES**

Plastics: 020 temp splint (0.5mm)

1.- Take an impression before prepare the teeth. Obtain the plaster model. You can work over all arch or over a special area.

2.- If a bridge will be mounted over a edentated area, sintetic teeth should be placed in the model to remplace the missing teeth. Can be fixed with wax. As in all techniques the model should be trimmed as deep as possible. A hole in the palatal area will help suction if necessary.

3.- Adapt the sheet

4.- Use a scissors or a knife to cut the material. And separate from model. Take care to avoid breaking the model while separation.

5.- If the temporary will be done by a laboratory, it should be sent to the clinc as it is now.

6.- Leave the adaptation to cover 1 or 2 teeth at both sides, and 35mm of soft tissue.

7.- While preparation it can be inserted in the mouth to get references of the borders. The adaptation is a transparent anatomical design that leaves a perfect occlusion.

8.- After trimmed, lubricate it. Place inside the powder and liquid to fill until the soft tissue height. If there are bubbles they can be seen easily now. When the resin starts to harden, place into the mouth.

9.- The exact quantity can be seen into the mouth. All must be covered without spaces or bubbles. When the resin starts to cure, take it out to avoid a bonding with the teeth.

10.- The curing temperature can be noted trough the plastic. It gets enough hardness to get out of mouth before the maximum curing temperature. This way we avoid discomfort to the patient.

11.- After taking out the mouth, leave to cure completely in heat water or under press.

12.- When the resin is completely cured, the plastic can be removed. The temporary will have shiny surface. Just need to trim the margins and widen the interdental spaces with a disc.

13.- Cement in the usual way

## **FLUORIDE TRAYS** **SURGERY APPLICATIONS**

Plastic 020 temp splint (0,5mm)

The adaptation is very thin well adapted to teeth and tissues. Does not alter occlusion neither movements. With this material can be obtained many other uses not said in this notes.

- 1.- The adaptation for customized fluoride trays is done in the usual way, and cut with scissors.
- 2.- Preparation of the model
  1. if the adaptation will be for retainer or to press any edentated area, we must take out first the teeth that will be extracted and with a knife we should prepare the area. Other way is to take a surgery impression, and get the model from the resulting. The space for the fluor or medicament is got by refractory mass before adapting.
  2. if the adaptation is to fix any medicin fixed into an area, as in periodontal procedures, some refractory mass should be placed over the model to leave the space.

NOTE: for best adaptation details, reduce the model height as much as possible, and make a hole in palatal area before adapting.

## **DUPLICATES**

Plastic: 020 copying

This technique can be used also for veneer crowns bases, for telescopic crowns, model duplications and duplicates of 1 only piece. As well as porcelain cores.

The copying material produces an optimal adaptation, dimensional stability and 100% calcination. The copies can be got without alterate the model. The core margins can be seen trough the plastic. The neck is fixed accurately due to the stability of the base. If windows are wanted on the veneers, they are easy got taking out the plastic.

The copying material is also suitable for temporary crown and bridges, because the autocure resin does not bonds to it.

- 1.- Place the cores over refractory mass as base. Use a small ball for each piece. If there are retentive areas, they should be covered with refrac. mass. We recommend to paint the cores with thick nail paint.
- 2.- Spray the cores mounted on the machine with silicone spray. Place the cores separated, and separated from the borders
- 3.- Adapt the plastic. In the heating process, the sheet goes up and then gets theyr initial height, and get a transparent aspect. This time you can adapt it. If you leave to drop, the adaptation will be too thin.
- 4.- After adapt, turn back the heater and leave the suction 30-60 sg. The plastic will get a translucent aspect while cools.
- 5.- If the retain areas have ben covered just pulling back, ca be separated from the plastic.
- 6.- Use a scissors to take out the material.
- 7.- The individual adaptations can be separated from the cores cutting the plastic with any heated instrument, no closer than 1mm under the neck. Before separating paint the margin with a pencil over the plastic to know the level to cut.
- 8.- The final cutting can be done with curved scissors

9.- Spray the core with silicone spray. Place the plastic over the core and with wax , finish it depending the final use desired as veneer, crown, duplicate, etc Any retention system is available.

10.- Use the regular procedure. Before fuse, spray with silicone to assure a smooth surface bubble free.

### **WHITENING TRAYS**

Plastics: 040 mouthguard (1.0mm), 060 Mouthguard (1.5mm), 080 mouthguard (2.0mm). If extra thin tray is desired 020 temp splint (0.5mm) or 020 copying can be used.

The special technique is indicated in the dentaflux whiteners as Bladen or Blanqueamiento dentatech

The technique is the usual. Be careful, the heated soft sheets bonds to everything. The common varnish is not indicated because reduces the plaster porosity and so the suction capacity.

Foam can be placed on the plastic when adapted (see custom trays, but use dry foam). This avoid rinsing the gel too soon, but it can be too thick and less estetic. And it is difficult to clean. Anyway the molars must be left without foam to assure fixation.

To leave space the teeth to be trated should be painted, for example with nail varnish to create space for the gel. LEAVING THE NECKS FREE.

The first cut must be done 5mm under necks so with the finishing leave at 2mm.

We recommend to keep the models to repeat the adaptations if lost

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