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## 6.2 Leaks? - Replace the o-rings

If the coupler leaks, replace the o-rings. Press the defective o-ring(s) to one side between two fingers until a loop is formed and then slide the rings to the front and remove them. Insert new o-ring(s) in the recesses and apply Lares One-Step conditioner.

 **Caution:** Do not use Vaseline or any other type of grease on O-rings. It may result in handpiece or coupler damage.

## 7.0 Care Instructions

### 7.1 Cleaning Exterior

The Lares *Fluoresce HD* LED Coupler can be cleaned on the outside with an isopropyl alcohol solution.

 **Caution:** Do not place the Lares *Fluoresce HD* LED Coupler in an ultrasonic cleaner since this may result in malfunction.

#### WARNING

##### **Hazard from non-sterile products.**

There is a risk of infection from contaminated handpieces and accessories.

- ▶ Take suitable personal protective measures (i.e. exam gloves, eye protection, 5 µm particulate filter mask).

- ▶ Remove all residual cement, composite or blood without delay.
- ▶ Clean and maintain the coupler as soon as possible after treatment.
- ▶ Remove the coupler from the dental unit hose.
- ▶ The coupler must be dry when transported for cleaning and sterilizing.
- ▶ Do not submerge or place it in any solution.

### 7.2 Cleaning

#### CAUTION

##### **Malfunctions from cleaning in the ultrasonic unit.**

Ultrasonic cleaning may cause severe damage to the product.

- ▶ Only clean manually.

#### 7.2.1 Manual Cleaning - external

Accessories required:

- Warm tap water 40 °C ± 5 °C (104 °F ± 10 °F)
- Brush, e.g. medium-hard toothbrush or comparable nylon-bristled dental instrument cleaning brush.

##### **Manual Cleaning Process**

Coupler disassembly not required.

- ▶ Rinse product under warm flowing tap water to remove gross oil.
- ▶ Brush off while rinsing under flowing warm tap water for 30 to 45 seconds to remove all soil. Brush in a direction parallel to the o-ring grooves and completely around the coupler.
- ▶ Use filtered pressurized air to completely dry coupler. If soil is still visible when dry, repeat all the cleaning steps.
- ▶ Apply Lares One-Step Handpiece Conditioner to the exterior of the coupler. Spray the conditioner over the coupler into a towel or sink.
- ▶ Remove excess conditioner from the exterior of coupler with a towel.

#### 7.2.2 Automated External Cleaning

Not applicable.

#### 7.2.3 Manual Cleaning of the Inside

Not applicable.

#### 7.2.4 Automated Internal Cleaning

Not applicable.

### 7.3 Disinfection

#### CAUTION

##### **Malfunctioning from using a disinfectant bath or disinfectant containing chlorine.**

Handpiece damage may result.

- ▶ Never disinfect in a thermo-disinfector.

#### 7.3.1 Manual Disinfection - External

The Lares *Fluoresce HD* LED Coupler can be disinfected on the outside with chemical disinfectants (spray disinfectant or wipe clean disinfectants). The coupler should not be immersed in any disinfecting solutions.

#### 7.3.2 Manual Disinfection - Internal

Not applicable.

#### 7.3.3 Machine Disinfection - External & Internal

Not applicable.

### 7.4 Drying

Manual Drying

- ▶ Blow off the outside and inside with compressed air until water drops are no longer visible.

#### 7.4.1 Automatic Drying

Not applicable

### 7.5 Care Products and Systems

#### - Servicing

#### CAUTION

##### **Premature wear and malfunctions from improper servicing and care.**

Reduced product life.

- ▶ Perform proper care regularly!

##### **Note**

Lares Research only guarantees that its products will function properly when the care products used are those listed as accessories, as they were tested for proper use on our products.

#### 7.5.1 Care Products & Systems - Servicing: Care with Lares One-Step Handpiece Conditioner

Lares Research recommends applying handpiece conditioner after each time it is used, i.e. after each cleaning and before each sterilization.

- ▶ Remove LED lamp assembly.

#### CAUTION

##### **Malfunction from improper servicing and care.**

Reduced product life or damaged device.

- ▶ Be sure to remove LED Lamp assembly prior to autoclaving.

#### 7.5.2 Inspection & Function Testing

*Fluoresce HD* LED MX Couplers: Visually inspect for damage and wear. Service or discard damaged couplers.

#### 7.5.3 Packaging Prior to Autoclaving

##### **Note**

The autoclave bag must be large enough for the coupler so that the bag is not stretched.

The quality and use of the autoclave bag must satisfy applicable standards, be cleared to market and be suitable for the autoclave procedure!

- ▶ Seal each coupler individually in an autoclave bag!

#### 7.5.4 Autoclaving

Process couplers in a steam sterilizer (autoclave) in accordance with EN 13060/ISO 17665-1.

#### CAUTION

##### **Premature wear and malfunctions from improper servicing and care.**

Reduced product life.

- ▶ Before each autoclave cycle, service the coupler with Lares One-Step Handpiece Conditioner.

#### CAUTION

##### **Contact corrosion due to moisture.**

Leaving coupler in autoclave after cycle has finished may damage the coupler.

- ▶ Immediately remove the product from the steam sterilizer after the sterilization cycle to reduce condensation!

#### CAUTION

##### **Do not use Flash Sterilization Techniques.**

These techniques are not validated for efficacy and damage to the coupler may result.

The Lares Research *Fluoresce HD* LED Coupler has a maximum temperature resistance up to 135 °C (275 °F).

Autoclave with gravitation process for at least 15 minutes at 132 °C ± 1 °C (270 °F ± 1.8 °F).

Drying time: 15 minutes minimum.

When autoclaving multiple couplers in one autoclave cycle ensure that the autoclave's maximum load is not exceeded. Follow the autoclave manufacturer's Instructions for Use.

#### 7.5.5 Handling & Storage After Autoclaving

- ▶ Keep couplers in the autoclave bag until next use.
- ▶ Cleaned and autoclaved couplers should be stored protected from dust with minimum exposure to germs in a dry, dark and cool space.
- ▶ Re-install LED Lamp assembly.

##### **Note**

The reprocessing instructions provided have been validated by Lares Research as being CAPABLE of preparing a *Fluoresce HD* LED Coupler for re-use. It remains the responsibility of the practitioner to ensure that the processing was actually performed using equipment, materials and personnel in the practitioner's facility to achieve the desired result. This requires validation and routine monitoring of the process as specified by the autoclave manufacturer. Use only FDA cleared to market equipment and materials for autoclaving and routine monitoring of the autoclaving process.

#### WARNING

##### **Hazard from cross infection.**

There is a risk of cross infection from patient to patient to dental staff from contaminated couplers.

- ▶ Do not deviate from the specified coupler reprocessing steps (i.e. cleaning, autoclaving and storage).

## 8.0 Accessories

### Replacement Parts

Description	Part ID
Wrench	10363
Gasket	10905
Replacement <i>Fluoresce HD</i>	11515
LED Lamp	
O-ring Kit	10074
One-Step Handpiece	10083
Conditioner	

## 9.0 Limited Warranty:

Each Lares *Fluoresce HD* LED Coupler (and LED lamp assembly) is warranted against defects in materials and workmanship for a period of 1 year from the date of purchase.

### 9.1 Additional Conditions of Warranty:

1. Warranty registration is automatic as of shipping date (Outside the US warranty registration may be required).
2. The *Fluoresce HD* LED Coupler must be operated and maintained in accordance with procedures outlined in these instructions.
3. The *Fluoresce HD* LED Coupler must not be subjected to abuse or neglect.
4. The *Fluoresce HD* LED Coupler must not be repaired or disassembled by anyone other than Lares Research or your authorized Lares distributor. Lares Research will repair or replace, at its discretion, without charge any defective parts covered by this warranty provided the swivel coupler is returned to the factory, transportation prepaid. (Outside the US return to your authorized Lares distributor.) Lares Research makes no other warranties expressed or implied.