MTA | BIOSEAL

**LIVRO
BRANCO**

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
Dados de literatura científica

- [1] Endodontic Practice – L. I. Grossman - 10th Edition - 1981
- [2] Physicochemical properties and surfaces Morphologies evaluation of MTA Fillapex and AH Plus – Borges *et al.* – The scientific world Journal – 2014
- [3] A mineral trioxide aggregate sealer stimulated mineralization – Gomes-Filho *et al.* – Journal of Endodontics – 2009
- [4] Microscopic appearance and apical seal of root canals filled with gutta-percha and ProRoot Endo Sealer after immersion in a phosphate-containing fluid – Weller *et al.* – International Endodontic Journal – 2008
- [5] Mineral trioxide aggregate: a comprehensive literature review-Part III: Clinical applications, drawbacks, and mechanism of action - Parirokh M, Torabinejad M – Journal of Endodontics - 2010
- [6] Internal documentation – ITENA CLINICAL
- [7] Investigation of Mineral Trioxide Aggregate for root end filling in dogs, Torabinejad et al. – Journal of Endodontics – 1995
- [8] Biototoxicity of commonly used root canal sealers: a meta-analysis – Kaur *et al.* – Journal of Conservative Dentistry – 2015
- [9] Ion Release, Porosity, Solubility, and Bioactivity of MTA Plus Tricalcium Silicate, Gandolfi *et al.* – Journal of Endodontics – 2014
- [10] Chemical characteristics of mineral trioxide aggregate and its hydration reaction – Chang – Restorative Dentistry and Endodontics - 2012
- [11] Evaluation of cytotoxicity and physicochemical properties of calcium silicate-based endodontic sealer MTA Fillapex - Silva *et al.* – Journal of Endodontics – 2013
- [12] Cellular mechanisms of bone remodelling – Erik Fink Eriksen – Reviews in Endocrine and Metabolic Disorders – 2010
- [13] EM evaluation of bacterial biofilm and microorganisms on the apical external root surface of human teeth – Leornado *et al.* – Journal of Endodontics – 2002
- [14] Pathogenesis of apical periodontitis and the causes of endodontic failures – Critical Reviews in Oral Biology and Medicine Nair PN – 2004
- [15] Dental biofilm: difficult therapeutic targets - Socransky and Haffajee – Periodontology 2000 – 2002
- [16] Antibiofilm activity, pH and solubility of endodontic sealers – Faria-Júnior *et al.* – International Endodontic Journal – 2013
- [17] Reaction of rat connective tissue to implanted dentin tubes filled with mineral trioxide aggregate or calcium hydroxide – Holland *et al.* – Journal of Endodontics – 1999
- [18] MTA and F-doped MTA cements used as sealers with warm gutta-percha. Long term study of sealing ability - Gandolfi and Prati – International Endodontic Journal – 2010

- [19] Push-out strength of modified Portland cements and resins – Ianoco *et al.* – American Journal of Dentistry – 2010
- [20] Effect of MTA-based sealer on the healing of periapical lesions – Gomes-Filho *et al.* – Journal of Applied Oral Science – 2013
- [21] pH changes in dental tissues after root canal fillin with calcium hydroxide – Tronstad *et al.* – Journal of Endodontics – 1981
- [22] Characterization of calcium oxide in root perforation sealer materials – Estrela *et al.* – Brazilian Dental Journal – 2012
- [23] Mineral trioxide aggregate (MTA) solubility and porosity with different water-to-powder ratios – Fridland and Rosado – Journal of Endodontics – 2003
- [24] *In Vitro* Biocompatibility, inflammatory response, and osteogenic potential of 4 root canal sealers: Sealapex, Sankin Apatite Root Sealer, MTA Fillapex and Iroot SP root canal sealer – Chang *et al.* – Journal of Endodontics – 2014
- [25] Odontoblast-like cytodifferentiation of human dental pulp cells in vitro in the presence of a calcium hydroxide-containing cement – Seux *et al.* – Archive of Oral Biology - 1991
- [26] Mineral Trioxide Aggregate-based endodontic seamer stimulates hydroxyapatite nucleation in human osteoblast-like cell culture – Salles *et al.* – Journal of Endodontics – 2012
- [27] Rat tissue reaction to MTA Fillapex – Gomes-Filho *et al.* – Dental Traumatology – 2011
- [28] Effects of calcium hydroxide-containing pulp-capping agents on pulp cell migration, proliferation, and differentiation – Schröder U. – Journal of Dental Research – 1985
- [29] Cytotoxic effects of four different root canal sealers on human osteoblasts – Jung *et al.*- PLOS One – 2018
- [30] Interfacial adaptation of adhesive materials to root canal dentin – Perdigao *et al.* – Journal of Endodontics – 2007
- [31] Depth and percentage of penetration of endodontic sealers into dentinal tubules after root canal obturation using a lateral compaction technique: a confocal laser scanning microscopy study – Ordinola-Zapata *et al.* – Oral Surgery, Oral Medicine, Oral Pathology, Oral Radiology and Endodontics – 2009
- [32] Filling effectiveness and dentinal penetration of endodontic sealers: a stereo and confocal laser scanning microscopy study – Silva RV *et al.* – Brazilian Dental Journal – 2015
- [33] The influence of the smear layer on dentinal tubule penetration depth by three different root canal sealers: an in vitro study – Kokkas *et al.* – Journal of Endodontics – 2004
- [34] Penetration of dentinal tubules by endodontic sealer cements in extracted teeth and in vivo – Mamootil and Messer – International Endodontic Journal – 2007
- [35] Influence of the calcium hydroxide intracanal dressing on dentinal tubule pentration of two root canal sealers – Cruz *et al.*- European Endodontic Journal – 2017

- [36] Comparative evaluation of dentinal penetration of three different endodontic sealers with and without smear layer removal – scanning electron microscopic study – Sonu *et al.* – Saudi Endodontic Journal – 2016
- [37] Radiopacity and flow of different endodontic sealers – Tanomaru-Filho *et al.*- Acta Odontologica Latinoamericana – 2013
- [38] Physical properties of 5 root canal sealers – Zhou *et al.* – Journal of Endodontics – 2013
- [39] Physical properties of MTA Fillapex sealer – Vitti *et al.*- Journal of Endodontics – 2013
- [40] Dissolution of a Mineral Trioxide Aggregate sealer in endodontic solvents compared to conventional sealers – Alzraikat *et al.* – Original Research Dental Materials – 2015
- [41] Evaluation of the pH, calcium release and antibacterial activity of MTA Fillapex – Kuga *et al.*- Revista de Odontologia Da Unesp – 2013
- [42] Removing root canal obturation materials: a comparison of rotary file systems and re-treatment agents – Ring *et al.* – Journal of the American dental association – 2009
- [43] Retreatability of 2 Mineral Trioxide Aggregate-based root canal sealers: a cone-beam computed tomography analysis – Neelakantan *et al.* – Journal of Endodontics – 2013