NUMERYS GF

WHITE PAPER

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1. Schwartz RS, Robbins JW. Post placement and restoration of endodontically treated teeth: A literature review. *J Endod.* 2004;30:289–301
2. Koutayas SO, Kern M. All-ceramic post and cores: The state of the art. *Quintessence Int.* 1999;30:383–92.
3. Aquaviva S, Fernandes and Coutinho. Factors determining post selection: A literature review. *J Prosthet Dent* 2003;90:556-62.
4. Sorrentino, R., Di Mauro, M.I., Ferrari, M. et al. Complications of endodontically treated teeth restored with fiber posts and single crowns or fixed dental prostheses-a systematic review. *Clin Oral Invest* (2016) 20: 1449.
5. Afroz S, Tripathi A, Chand P, Shanker R. Stress pattern generated by different post and core material combinations: A photoelastic study. *Indian J Dent Res* 2013;24:93-7.
6. M. Ona, N. Wakabayashi, T. Yamazaki, A. Takaichi, and Y. Igarashi, "The influence of elastic modulus mismatch between tooth and post and core restorations on root fracture," *International Endodontic Journal*, vol. 46, no. 1, pp. 47–52, 2013
7. Jianliang PANG et al. Fracture behaviors of maxillary central incisors with flared root canals restored with CAD/CAM integrated glass fiber post-and-core. *Dent Mater J* 2019; 38(1): 114–119
8. Aashwini Lamichhane, Chun Xu1, Fu-qiang Zhang. Dental fiber-post resin base material: a review. *J Adv Prosthodont* 2014;6:60-5
9. Boudrias P, Sakkal S, Petrova Y. Anatomical post design meets quartz fiber technology: rationale and case report. *Compend Contin Educ Dent.* 2001 Apr;22(4):337-40, 342, 344.
10. Dyer et al. Effect of fiber position and orientation on fracture load of fiber-reinforced composite. *Dental Materials* (2004) 20, 947–955
11. Zicari et al. Mechanical properties and micro-morphology of fiber posts. *Dental Materials* 29 (2013) e45-e52.
12. Perdigão J, Gomes G, Lee IK. The effect of silane on the bond strengths of fiber posts. *Dent Mater* 2006;22:752-8
13. Ferrari et al. Long-term retrospective study of the clinical performance of fiber posts. *Am J Dent*, 2007 Oct;20(5):287-91
14. Kremeier K, Fasen L, Klaiber B, Hofmann N. Influence of endodontic post type (glass fiber, quartz fiber or gold) and luting material on push-out bond strength to dentin in vitro. *Dent Mater.* 2008 May;24(5):660-6.
15. Muttlib et al. Intracanal Adaptation of a Fiber Reinforced Post System as Compared to a Cast Post-and-Core. *Acta stomatol Croat.* 2016;50(4):329-336.
16. Anonymous. The glossary of prosthodontic terms. *J Prosthet Dent.* 2005 Jul;94(1):10-92
17. Turner CH. Post-retained crown failure: a survey. *Dent Update.* 1982 May;9(4):221, 224-6, 228-9 passim
18. Musikant BL, Deutsch AS. Endodontic posts and cores. Part II. Design of the Flexi-post. *J Ala Dent Assoc.* 1985 Fall;69(4):42-6
19. Geramipannah F, Rezaei SMM, Sichani SF, Sichani BF, Sadighpour L. Microleakage of Different Post Systems and a Custom Adapted Fiber Post. *J Dent (Tehran).* 2013 Jan;10(1):94-102
20. Balkenhol et al. Survival time of cast post and cores: A 10-year retrospective study. *J Dent.* 2007 Jan;35(1):50-8. Epub 2006 Jun 5.
21. Fraga RC, Chaves BT, Mello GS, Siqueira JF, Jr. Fracture resistance of endodontically treated roots after restoration. *J Oral Rehabil.* 1998 Nov;25(11):809-13.
22. Christensen GJ. Why all ceramic crowns? *J Am Dent Assoc.* 1997;128:1433–5
23. Blatz MB. Long term clinical success of all-ceramic posterior restorations. *Quintessence Int.* 2002;33:415–26
24. Liu et al. Use of a CAD/CAM-fabricated glass fiber post and core to restore fractured anterior teeth: A clinical report. *J Prosthet Dent.* 2010 Jun;103(6):330-3
25. ITENA Clinical – Internal documentation – Bench test - 2018

26. Ruschel et al. Mechanical properties and superficial characterization of a milled CAD-CAM glass fiber post. *J Mech Behav Biomed Mater.* 2018 Jun;82:187-192
27. ITENA Clinical – NUMERYYS GF technical drawings - 2018
28. ITENA Clinical – NUMERYYS GF Radiopacity assessment - 2018
29. ITENA Clinical – DENTOCLIC glass fiber Technical Paper – 2018
30. Garbin et al. Biomechanical behaviour of a fractured maxillary incisor restored with direct composite resin only or with different post systems. *International Endodontic Journal*, 43, 1098–1107, 2010
31. Spazzin et al. Influence of Post and Resin Cement on Stress Distribution of Maxillary Central Incisors Restored with Direct Resin Composite. *Operative Dentistry*, 2009, 34-2, 223-229
32. ITENA Clinical – Test Report: Mechanical Performances -2018
33. Vano et al. The adhesion between fibre posts and composite resin cores: the evaluation of microtensile bond strength following various surface chemical treatments to posts. *International Endodontic Journal*, 39, 31–39, 2006
34. Galbano et al. Evaluation of the Flexural Strength of Carbon Fiber-, Quartz Fiber-, and Glass Fiber-Based Posts. *JOE — Volume 31, Number 3, March 2005*
35. K.-J.M. Söderholm and S.-W. Shang. Molecular Orientation of Silane at the Surface of Colloidal Silica. *J Dent Res* 72(6):1050-1054, June, 1993