Protocolo de CAD

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The screenshot displays the worknc Dental software interface. At the top, there is a navigation bar with icons for various dental procedures: ACCUEIL, PRÉPARATION, POST-TRAITEMENT, IMPORT, GAMME, CAVITE, IMPLANTS, ORIENTATION, ATTACHE, REPRISE, and USAGE. The 'IMPORT' icon is highlighted. The main workspace shows a 3D model of a dental implant assembly within a transparent wireframe bounding box. On the left, there is a sidebar with an 'Import' section containing a file named '1 - 82799_20190211_1025_TRIOS_Basic (TRIO...'. Below this, there are dropdown menus for 'Type', 'Client', and 'Tag personnalisé', all currently empty. A 'SUIVANT' button is located at the bottom left. On the right side, there is a vertical toolbar with various manipulation icons, including a 'CAD' icon. The 'worknc Dental' logo is visible in the top right corner.

The screenshot displays the worknc Dental software interface. At the top, a navigation bar includes icons for 'ACCUEIL', 'PRÉPARATION', and 'POST-TRAITEMENT'. A secondary bar contains icons for 'IMPORT', 'GAMME', 'CAVITÉ', 'IMPLANTS', 'GÉNÉRATION', 'ATTACHE', 'REPRISE', and 'USAGE'. The 'GAMME' icon is highlighted. The main workspace shows a 3D model of a dental crown on a stem, rendered in green, positioned within a wireframe bounding box. A grey semi-transparent model of the crown's internal structure is visible behind it. On the right side, a vertical toolbar contains various icons for navigation and editing. On the left, a sidebar panel shows a list of items under 'Gamme', with 'LYRA_T10_V19_open' selected. Below this, technical specifications are listed: 'Qualité de finition : Moyenne', 'Conçue pour : Couronne anatomique', 'Nécessite : Ligne marginale requise', and 'Outils: T20 Ball D1.13;'. A 'SUIVANT' button is located at the bottom left of the sidebar.

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The screenshot displays the worknc Dental software interface. At the top, a green navigation bar contains icons for various stages: IMPORT, GAMME, CAVITÉ, IMPLANTS (highlighted), ORIENTATION, ATTACHE, REPRISE, and USAGE. The 'worknc Dental' logo is on the right. On the left, a sidebar shows the file name '1 - 82799_20190211_1025_TRIOS_Basic_(TRIO...' and a 'SUIVANT' button at the bottom. The main workspace shows a 3D model of a dental implant with a grey abutment and a white crown, positioned within a blue wireframe box. A vertical toolbar on the right contains various manipulation tools. Below the main workspace, a panel labeled 'Affichage' shows icons for different views, and a text prompt reads 'Sélectionner une entité pour la modifier'.

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