# PURE OFFICE

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## EN DIRECTIONS FOR USE

**PURE OFFICE** is a bleaching gel intended for medical in-office use only. The material is used by the dentist to whiten individual teeth. **PURE OFFICE** is ideally suited to whiten devital teeth including endodontically treated teeth (intracanal use). The 4:1 automix syringe contains bleaching gel and activator. The dual syringe contains 5 g of material. The hydrogen peroxide content of the ready-mixed gel is 35%.

### INDICATIONS

- Intracanal whitening of discolored, devital teeth
- Intracanal whitening of discolored teeth caused by disease, by injury or iatrogenically, as an alternative to a crown or veneer
- Intracanal treatment of tooth discolorations caused by prescription drugs (e. g. tetracycline)

### CONTRAINDICATIONS

- Exposed dentin
- Defective filling margins
- Hypersensitive patients
- Allergy to or intolerance of the ingredients
- Pregnancy and lactation (due to insufficient data)

Some patients may experience slight sensitization in their teeth during or shortly after treatment, but this normally subsides after a short time.Stop treatment immediately if undue pain is being experienced.

### PREPARATIONS

Determine the initial tooth color. Use air flow or a cleansing paste to clean the patient's teeth. Use a cheek / lip retractor or similar instrument. Protect the gingiva (e. g. with PURE PROTECT or by applying a rubber dam). Do not anaesthetize the patient. The separate components (gel and activator) as well as the mixed bleaching gel irritate the skin and the mucous membranes, and cause severe damage to the eyes. Wear protective gloves/ protective clothing/protective goggles/ face mask. IN CASE OF EYE CONTACT: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a doctor.

### APPLICATION

First remove both the fixing clamp and the sealing cap, facing away from the patient and from yourself. Replace the sealing cap with a mixing tip. Now squeeze the gel from the syringe until a homogeneously colored mass with gel-like consistency is coming out of the tip. Discard the initial inhomogeneous mass. Storage for long periods at high temperatures may sometimes cause bubbles to form in the gel. This can also lead to inhomogeneous and more liquid mixtures, especially in the last quarter of the dual syringe. If this occurs, do not use the affected parts of the material. The gel should be applied in a thickness of about 1- 2 mm to the labial surfaces of the teeth to be bleached. The point of the mixing tip is designed to spread the gel evenly. To enhance the effect, you can use a standard plasma arc light or laser. In this case, both the dentist and patient should wear safety goggles. Please follow the manufacturer's instructions. Do not overheat. After a reaction time of about 10 minutes, aspirate the gel and wipe off the residue. If the bleaching effect is insufficient, the procedure can be repeated twice until the desired effect has been achieved. Before each treatment, ensure that the gingiva protection is still sufficient. Then determine the tooth color again. After the final bleaching procedure, aspirate the gel, wipe it off and then rinse thoroughly with water. Now dry the tooth surface and polish with polishing paste. Subsequently, a remineralization with PURE CARE desensitizing gel or with a fluoride solution or fluoride gel is recommended. This is also suitable for treating pain if sensitization occurs. If there is still enough material for another bleaching treatment, remove the mixing tip and replace the blue cap. Store the material in a refrigerator. Use a new mixing tip for the next treatment. The patient should not eat, drink or smoke for 1 hour after treatment. Substances with strong colors (e. g. red wine, coffee) and foods containing acid should also be avoided in the days immediately following the treatment.

Note: Any composite fillings and restorative materials made of ceramic or metal are not bleached. Fillings should there fore be replaced 2 weeks after treatment.

### INTRACANAL TREATMENT OF NON-VITAL TOOTH

The gel is also applied through an access opening on the lingual side of the tooth. First, check the height of the alveolar ridge. Then completely remove the filling material from the pulp chamber.

Only apply one layer of glass ionomer filling with a thickness of 1- 2 mm to the bottom of the pulp chamber on the same level with the alveolar ridge or coronally. Apply the gingiva protection to both the labial and the lingual sides of the tooth. Now, apply the bleaching gel to the labial surface as described under „Application” above. Additionally, apply the gel through the access opening. Let the gel act for 10 min. Afterwards, carefully aspirate the gel and lightly rinse with water. If the bleaching effect is insufficient, this procedure can also be repeated twice until the desired effect has been achieved. After the final bleaching procedure, aspirate the gel and rinse thoroughly with water. Carefully remove the gingiva protection. Next, close the access opening with a temporary filling material which does not contain eugenol. Wait two weeks before preparing the permanent filling.

### SHELF LIFE / STORAGE

Do not use PURE OFFICE after the expiry date, because an adequate bleaching effect is no longer guaranteed. Refrigerate the bleaching gel at 3 - 9 °C. Do not freeze.The product contains hydrogen peroxide, potassium nitrate, sodium fluoride, sodium hydroxide, glycol derivatives, thickening agent and coloring.

### TROUBLESHOOTING

Insufficient bleaching effect	Repeat the bleaching procedure (several times if necessary)
Patient is experiencing pain	The gel has come into contact with exposed dentine or a loose filling (please cover before treatment). The patient is particularly sensitive (depending on the individual tooth structure).
There are white spots on the enamel after the bleaching procedure	Local demineralization or dehydration. Reversible. Disappears after fluoridation.
Gingiva is colored white after bleaching procedure can be reversed)	The gingiva is not adequately protected, or the protection is loose
The gel from the cartridge has an in- homogeneous color (colorless / red)	Air bubbles in the gel as a result of being stored at too high a temperature or for too long.
Reduced composite-enamel adhesion after bleaching procedure	Peroxide residue. Leave for 1- 2 weeks before restoration.

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