



# DENTO-VIRACTIS INSTRUMENTS 55

Dilution:  
0,5%



<b>Description</b>	Dento-Viractis 55 is recommended for <b>cleaning and pre-disinfection by immersion of dental instruments</b> , after use and before sterilisation or disinfection.
<b>Scope of application</b>	For use on heat-resistant and heat-sensitive reusable medical instruments and devices.
<b>Characteristics</b>	<p><b>Prevents biofilm formation</b> and removes dirt</p> <ul style="list-style-type: none"> <li>✓ Transforms hard water into soft water with its powerful sequestering agent: <b>your instruments are revived!</b></li> <li>✓ <b>Enzyme preserved</b> thanks to the powdered product presentation.</li> <li>✓ <b>Reduces the risk of cross-contamination:</b> the rejected bath is not con contaminated</li> </ul>
<b>Commercial presentation</b>	2 kg bucket with measuring scoop 5 kg bucket with measuring scoop Bucket of 50 water-soluble doses of 15 g

Antimicrobial activity			Concentration	Contact time
<b>Bactericidal</b>	EN 1040	Pseudomonas aeruginosa, Staphylococcus aureus	0,1%	15 min
	EN 13727	Pseudomonas Aeruginosa, Staphylococcus aureus, Enterococcus hirae, Acinetobacter baumannii	0,5%	
<b>Yeasticidal</b>	EN 1275	Candida albicans	0,5%	15 min
	EN 13624	Candida albicans	0,25%	
<b>Virucidal</b>	EN 14476	Herpes virus	0,5%	1 min
		Hepatite B		15 min
		Hepatite C		
		HIV -1		

DENTO-VIRACTIS



# DENTO-VIRACTIS INSTRUMENTS 55

# DENTO-VIRACTIS

<b>Composition</b>	Cationic and non ionic surfactants, Sequestering agents (phosphates), Enzyme Proteolytic, Perfume.
<b>Instructions for use</b>	Keep the container tightly closed and protected from moisture, in its original packaging and labelling. The 0.5% diluted solution is non-irritating to skin eyes. Do not use on aluminium. PH of the 0.5% solution in hard water: 10.3. <b>Preparation of the 0.5% solution:</b> <ol style="list-style-type: none"><li>1. Fill a graduated dip tank with a base of mains water at 30- 40°C.</li><li>2. Pour in the quantity of powder necessary to obtain the ready-to-use solution 5g/L (e.g. 25g for 5L). Close the bucket.</li><li>3. Stir until completely dissolved.</li><li>4. Make up to the desired volume with tap water at room temperature (20°C) until the desired volume is reached.</li></ol> <b>Use:</b> <ol style="list-style-type: none"><li>1. Completely immerse open or dismantled instruments and medical medical devices, taking care to irrigate the channels.</li><li>2. Contact time: 15 minutes.</li><li>3. Rinse thoroughly with mains water.</li><li>4. Renew the preparation as soon as it is visibly soiled or at the latest every 24 hours.</li></ol>
<b>Precautions for use</b>	Before use, read the label and product information. <b>DANGER.</b> ETHOXYL ALCOHOL C9 C11. Causes skin irritation. Causes severe eye damage. Harmful to aquatic organisms, causes long-term adverse effects. Avoid release to the environment. Wear protective gloves/protective clothing/eye/face protection. <b>IN CASE OF CONTACT WITH THE SKIN:</b> Wash with plenty of water. <b>IN CASE OF CONTACT WITH THE EYES:</b> Rinse cautiously with water for several minutes. Remove contact lenses if worn and if they can be easily removed. Continue rinsing. Call a POISON CENTER or doctor immediately. Dispose of contents and container as hazardous waste. CE mark obtained in 2019
<b>Labelling of DV55 concentrate</b>	

CE0459 mark obtained in 2019. Class IIb DM.  
Professional use  
Data sheet version: 1.0 of 09/2021

Distribué par :  
**CADENCE**  
2 bis, chemin du loup  
93290 Tremblay-en-France  
FRANCE



**SODEL - ISO 9001/13485**  
190 rue René Barthélemy, Pôle d'activités  
de l'Espérance, 14100 Lisieux, FRANCE.  
Tel: +33 (0)2 31 31 10 50  
[www.sodel-sa.eu](http://www.sodel-sa.eu) - [www.alkapharm.fr](http://www.alkapharm.fr)