




DENTO-VIRACTIS INSTRUMENTS 53

DENTO-VIRACTIS

Description	<p>Dento-Viractis 53 is recommended for the immersion disinfection of dental dental instruments, directly after pre-disinfectant and deproteinizing cleaning with DV55.</p> <div style="text-align: right;">  <p>Dilution: 2% ou 5%</p> </div>
Fields of application	<p>For use on heat-sensitive reusable medical instruments and devices.</p>
Characteristics	<p>Available in concentrate (1 L) and ready-to-use solution (5 L).</p>
Commercial presentation	<p>The bottle of 1 liter concentrated with graduated doser The can of 5 liters of solution ready to use</p>

Antimicrobial activity. Cleanliness condition - 20°C		Concentration	Contact time
Bactericide	EN 13727	Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae	0,25%
	EN 14561	Pseudomonas Aeruginosa, Staphylococcus aureus, Enterococcus hirae, Acinetobacter baumannii	1%
	IRM	MRSA, IRPA	2%
Fungicide	EN 13624	Candida albicans	0,1%
	EN 14562	Candida albicans	2%
	GERMANDE	Penicillium verrucosum	
	IRM	Aspergillus brasiliensis (niger)	2%
Virucide	EN 14476	Rotavirus, Herpes virus	0,5%
		Poliovirus	2%
		Adénovirus, Norovirus, HBV, HBC	2%
	EN 17111	Norovirus, Adénovirus, Vaccine**	1%
	IPL	HIV	0,5%
Mycobactericide	IRM	Mycobacterium terrae	2%
	CHEMILA (EN 14563*)	Mycobacterium terrae, Mycobacterium avium	5%



DENTO-VIRACTIS INSTRUMENTS 53

DENTO-VIRACTIS

Sporicidia	NF T 72-230	Bacillus subtilis	5%	60 min
	IRM	Bacillus subtilis	2%	15 min
	GERMANDE			
	CHEMILA (EN 14347***)	Bacillus subtilis, Bacillus cereus	5%	

The MRI, Germande and Chemila studies use **DV55** 0.5%/15min before **DV53**.
DV53 PAE corresponds to a 5% dilution of DV53 concentrate.

*According to the protocol adapted from EN 14563

** The Vaccine Virus is the representative virus of the enveloped viruses presented in the annexes of EN 14476+A2 (2019), EN 16777 and EN 17111, such as: Hepatitis B virus (HBV); Hepatitis C virus (HCV); Human immunodeficiency virus (HIV); Herpesviridae (Herpes virus); Coronavirus; Influenza, rabies, rubella, measles viruses;

*** According to protocol adapted from EN 14347

Composition	Disinfectant : Pentane-1,5-dial purified, Non-ionic and cationic surfactants cationic surfactants, Perfume.
Instructions for use	<p>Keep the container tightly closed in a dry and well-ventilated place, in its original packaging and labeling.</p> <ul style="list-style-type: none"> • Measuring bottle: squeeze the bottle to obtain a dose (from 2%). • 5L ready-to-use solution: pour the solution into a disinfection container with a drip tray and a lid. <p>Prepare the solution:</p> <ol style="list-style-type: none"> 1. Use a graduated disinfection tray with a lid. 2. Concentrate: Squeeze the sides of the bottle to fill the dosing chamber. Pour the required dose into the tray: 20mL/L (2% solution) for low and intermediate level disinfection, 50mL/L (5% solution) for high level disinfection. 3. Top up with water at room temperature (20°C) until the desired volume is reached. 4. The ready-to-use solution can be stored for up to 30 days in a closed container, protected from heat. It is recommended to renew it as soon as it is soiled or at the latest 7 days after its first use. <p>Usage:</p> <ol style="list-style-type: none"> 1. After pre-disinfectant cleaning with Dento-Viractis 55, check that medical devices are properly rinsed and drained. 2. Immerse open or dismantled instruments and medical devices completely, taking care to irrigate the channels for 15 minutes. 3. Rinse thoroughly, internally and externally, with a quality of water appropriate to the nature of the procedures for which the medical device is intended. 4. Dry with a drape (clean or sterile) and/or by insufflation of medical air.



DENTO-VIRACTIS INSTRUMENTS 53

DENTO-VIRACTIS

Precautions

DANGER - contains GLUTARAL. DIDECYL DIMETHYL AMMONIUM CHLORIDE.

Harmful if swallowed. Causes skin irritation. May cause skin allergy. Causes severe eye damage. Toxic by inhalation. May cause allergic or asthmatic symptoms or breathing difficulties by inhalation. May irritate the respiratory tract.

Very toxic to aquatic organisms, causes long-term adverse effects. Corrosive to the respiratory tract.

Avoid breathing vapors. Do not eat, drink or smoke when handling this product. Avoid release to the environment. Wear protective gloves, protective clothing, eye/face protection. Wear respiratory protective equipment.

IF ON SKIN: Wash with plenty of soap and water.

Wash with soap and water.

IF INHALED: Move person to fresh air and hold in a position comfortable for breathing.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easily removable. Continue rinsing.

IF IN EYES, IF SWALLOWED with discomfort, if skin irritation or rash occurs, if respiratory symptoms occur: call a POISON CENTER or doctor immediately.

Collect spilled material.

Dispose of contents and container as hazardous waste.

Do not discharge residue into drains or waterways.

Professional use.

Labeling of DV53 concentrate



Labeling DV53 PAE



CE0459 mark obtained in 2019. Class IIb DM. Data sheet version: 2.0 of 09/2021

Distribué par :

CADENCE

2 bis, chemin du loup
93290 Tremblay-en-France
FRANCE



SODEL - ISO 9001/13485

190 rue René Barthélemy, Pôle d'activités
de l'Espérance, 14100 Lisieux, FRANCE.
Tel: +33 (0)2 31 31 10 50
www.sodel-sa.eu - www.alkapharm.fr