




DENTO-VIRACTIS 57

ROTARY INSTRUMENTATION

DENTO-VIRACTIS

Ready to use

Description	<p>Dento- Viractis 57 is a product for the cleaning and disinfection of burs and rotary instrumentation prior to sterilisation.</p> <div style="text-align: center;">  </div>
Scope of application	<p>For use on contaminated reusable burs and rotary instruments.</p>
Characteristics	<ul style="list-style-type: none"> ✓ Dropout protein and mineral stains ✓ Contains an anti-corrosion agent ✓ Product ready for use ✓ Compatible with ultrasonic baths
Commercial Presentation	<p>The concentrated 2 liters bottle ready to use The concentrated 5 liters bottle ready to use</p>

Antimicrobial activity			Contact time
Bactericidal	EN 1040	Pseudomonas aeruginosa, Staphylococcus aureus	5 min
	EN 13727	Pseudomonas aeruginosa, Staphylococcus aureus,	15 min
	EN 14561	Enterococcus hirae	15 min
Fongicidal	EN 1275	Candida albicans,	15 min
	EN 13624	Candida albicans, Aspergillus Niger	15 min
	EN 14562		15 min
Virucidal	EN 14476	Herpès simplex virus	15 min
	EN 14476	Hepatite C virus (BVDV)	15 min
	Pasteur institute	Hepatite B virus	10 min



DENTO-VIRACTIS

ROTARY INSTRUMENTATION **57**

DENTO-VIRACTIS

Composition	Nonionic surfactant, disinfectant, perfume (allergenic fragrances: d-limonene and linalool).
Precautions	DANGER. POTASSIUM HYDROXIDE. ETHOXYLATED ALCOHOL C10. Flammable liquid and vapour. Toxic if swallowed. Causes severe skin burns and eye damage. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Keep container tightly closed. Do not breathe vapours. Wear protective gloves, protective clothing, eye protection, face protection. IF SWALLOWED : Immediately call a POISON CENTER or doctor. Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair) : Take off immediately all contaminated clothing. Rinse skin with water or shower. IF IN EYES : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Store in a well-ventilated place. Keep cool. Dispose of contents/container as hazardous waste. CE Marking obtained in 2019.
Instructions for use	Store container tightly closed, in its original packaging and labelling. Do not use on aluminum parts. Ready to use solution. 1. Totally immerse rotary instruments. 2. Contact time : 15 minutes. 3. After cleaning instruments, rinse with tap water 4. Renew the solution when visibly soiled or at least every 24 hours.

Labelling:



CE0459 Marking obtained in 2019. Class IIb DM. Version of the sheet: 2.0 of 09/2021

Distribué par :
CADENCE
2 bis, chemin du loup
93290 Tremblay-en-France
FRANCE



SODEL - ISO 9001/13485
190 rue René Barthélemy, Pôle d'activités
de l'Espérance, 14100 Lisieux, FRANCE.
Tel: +33 (0)2 31 31 10 50
www.sodel-sa.eu - www.alkapharm.fr