




# DENTO-VIRACTIS INSTRUGERM 54

Dilution:  
0,5%

<b>Description</b>	<p>Dento-Viractis Instrugerm 54 is recommended for the <b>pre-disinfection of medical devices.</b></p> 
<b>Scope of application</b>	Use on medical devices.
<b>Commercial presentation</b>	The 1 liter bottle with measuring device. The 5 liter can.

DENTO-VIRACTIS

		Antimicrobial activity	Concentration	Contact time
<b>Bactericidal</b>	EN 13727	Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae	0,5%	5 min
	EN 14561	Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae	0,5%	10 min
<b>Yeasticidal</b>	EN 13624	Candida albicans	0,5%	5 min
	EN 14562	Candida albicans	0,5%	10 min
<b>Virucidal</b>	EN 14476	Herpes virus	0,5%	5 min
		Hepatitis B		
		Hepatitis C	0,5%	10 min
	HIV			
EN 17111	Vaccinia virus	0,5%	10 min	



# DENTO-VIRACTIS INSTRUGERM 54

# DENTO-VIRACTIS

<b>Instructions for use</b>	Concentration: 0,5 % (5 ml/L) = (25ml/5L). Contact time: 10 minutes. Temperature: 20°C.
<b>Precautions for use</b>	See Material Safety Data Sheet. <b>DANGER.</b> CASE 94667-33-1: N, N-DIDECYL-N-METHYL-POLY (OXYETHYL) AMMONIUM PROPIONATE. Causes severe skin burns and eye damage. Causes severe skin burns and eye damage. Very toxic to aquatic organisms, causes long-term adverse effects. Long-term adverse effects. Wear protective gloves/protective clothing/eye protection. Wear protective gloves/protective clothing/eye protection/face protection. <b>IF ON SKIN (or hair):</b> Take off immediately all contaminated clothing. Remove all contaminated clothing immediately. Rinse skin with water/shower. Rinse skin with water/shower. <b>IF IN EYES:</b> Rinse cautiously with water for several minutes. Remove contact lenses if worn and if they can be easily removed. Continue rinsing. Call a POISON CENTER/physician immediately.
<b>Notice Instruction Manual (2020/08/03)</b>	1st marketing: 2020-07. Intended for exclusive use by professionals for the pre-disinfection of invasive medical devices, such as surgical instruments surgical instruments, operating room equipment, endoscopes, ophthalmology and contactology ophthalmology and contactology equipment. Do not use the product after its expiration date. Close the canister tightly after use. Do not reuse the pre-disinfection bath (or follow local local recommendations. Do not mix with other products. The packaging must be disposed of as hazardous waste as hazardous waste under the sole responsibility of the holder of the waste. of this waste. Do not dispose of residues in sewers or waterways. Any serious incident involving the device should be reported to the to the manufacturer and the competent authority. Check the compatibility of the treated material with the materials tested in the corrosion study COR010435_mat: compatibility with stainless steel, aluminum, HDPE, silicone, endoscope sheaths, tungsten, viton, brass and PVC.
<b>Labelling</b>	

CE0459 mark obtained in 2019. Version 1.0. of 09/2021

Distribué par :  
**CADENCE**  
2 bis, chemin du loup  
93290 Tremblay-en-France  
FRANCE



**FRANK LAB**  
3 avenue des Frênes  
Z.A. de l'Observatoire  
78180 Montigny-le-Bretonneux  
FRANCE