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fr

INSTRUCTION FOR USE

en

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es

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pt

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it

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de

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nl

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SOMMAIRE

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ULTRASONIC CLEANER

SUMMARY

INTRODUCTION	2
1 General notes on delivery	2
2 General instructions	2
3 Guarantee information	3
CHAPTER 1	4
1.1 Unit and packaging dimensions	4
CHAPTER 2	4
2.1 Installation	4
2.2 Safety	5
2.3 Electrical connection	5
CHAPTER 3	6
3.1 Description	6
3.2 Technical specifications	6
3.3 Ultrasound tank accessories	6
3.4 Advantages of ultrasound cleaning compared to manual cleaning	7
3.5 How ultrasound works	7
3.6 Ultrasonic cleaning methods	8
3.7 Hints on using the unit	8
CHAPTER 4	9
4.1 Instructions for use	9
4.2 Definition of the colours front panel	9
CHAPTER 5	9
5.1 Routine maintenance and precautions	9
5.2 Trouble shooting	10
CHAPTER 6	10
6.1 Scrapping instructions	10

All instruction in this INSTRUCTION FOR USE must be completely read before installation or operation.

INTRODUCTION

1. GENERAL NOTES ON DELIVERY

When the unit is delivered, check the state of the packaging (keep the packaging for future use). Open the pack and check the following:

- Check that the delivered unit corresponds to technical specifications;
- Check that there is no evident damage.

In the event of damage or missing parts, inform the shipping agent or dental equipment depot immediately, giving all details. None of the dimensions given in the manual are binding.

The drawings and all other documentation supplied with the unit are the property of the manufacturer. All rights reserved; distribution to third parties strictly prohibited.

The partial or total reproduction of any text or illustrations is strictly forbidden.

2. GENERAL INSTRUCTIONS

The purpose of this manual is to give instructions on:

- The installation of the unit;
- It's safe and efficient use;
- Routine maintenance.

The unit/equipment must be used in compliance with the instructions in this manual and exclusively for its intended use. The user is legally responsible for the installation, operation and use of the unit.

The manufacturer shall not be held responsible for any running problems or any breakages or injuries caused by the unit if it is used incorrectly or not maintained properly.

Intended use: this manual constitutes an integral part of the product and must be kept in an accessible place near the unit where it can be consulted easily and quickly. The unit is designed for use by professional, qualified persons only. The unit must be used for its intended purpose only

ULTRASONIC CLEANER: equipment designed and developed for the cleaning of hard-to-access surfaces: diamond cutters, endodontic files, hand tools, scissors, dental prostheses (crowns, inlays), orthodontic appliances and laboratory instruments. The unit must be used exclusively with acid or alkaline based aqueous cleaning solutions

In this manual symbols are placed beside certain descriptions, notes, etc. These symbols are used to attract the attention of readers to a particular note or explanation. Their meaning is explained below.

DESCRIPTION	SYMBOL
IMPORTANT SAFETY INFORMATION This symbol is used to draw the reader's attention to particularly important notions for operator safety.	
INFORMATION AND PRECAUTIONS This symbol refers to general indications and advice.	
STRICTLY FORBIDDEN This symbol means it is strictly forbidden to perform the operation in question. Non-observance may cause serious harm to the operator or damage to the equipment.	

3. GUARANTEE INFORMATION

The manufacturer guarantees the quality of its equipment, if used in accordance with the instructions supplied in this manual, and according to the conditions printed on the GUARANTEE CERTIFICATE.


Warning: the CUSTOMER USER must fill in all parts of the DETACHABLE SLIP of the guarantee certificate and send it.

The guarantee is effective from the delivery date of the unit to the customer, which is confirmed by the return of the guarantee slip correctly filled in and signed. In the case of dispute, the date indicated on the purchase invoice, which shows the serial number of the unit, will be considered valid.



Note: Please retain the original packaging and use it each time the unit needs to be transported. The use of different packaging could damage the product during transportation.

CHAPTER 1


1.1 Unit and packaging dimensions

1. Dimensions of tank with cover

2. Tank casing dimensions

Capacity in liters: 3L

CHAPTER 2

2.1 Installation

The ultrasound tank has been tested and approved in our factory.

Remove the packaging and install the unit, taking the following precautions:

- 1 Position the ultrasound tank on a flat surface and respect a clearance of at least 50-80 mm on all sides.
- 2 Do not install the tank close to a sinks, washbasins, etc., to avoid contact with water and possible electrical short circuits.
- 3 Do not install close to any heat source.
- 4 Install in a normally ventilated room.
- 5 Make sure there are no kinks in the power supply lead and that there are no obstacles between the lead and the electricity socket.
- 6 Plug in to the socket; the unit is ready for use.

DESCRIPTION	SYMBOL
In order to safeguard the operator, the unit and the surrounding area, remember, when choosing the place of installation of any ultrasound equipment, to take into account the effects of the gradual evaporation and vaporisation of the cleaning fluid, especially if the unit is used without the cover. ALWAYS USE THE UNIT WITH THE COVER ON. DO NOT INSTALL THE APPLIANCE IN THE "PATIENT AREA".	
DANGER OF CONTAMINATION After use, the cleaning solution may contain contaminating residues: wear specific PPE when washing and cleaning the appliance and replacing the cleaning solution. The responsible authority must teach the operator how to use the unit safely.	
If particularly volatile and aggressive solutions are used, the fumes and vapours may spread into the surrounding environment and cause problems for any equipment and structures situated nearby.	
Do not use the unit for disinfecting, cleaning or rinsing contact lenses.	

2.2 Safety

Read this safety information carefully before using. The nonobservance of this information may cause accidents or damage to the equipment.

DO NOT USE INFLAMMABLE LIQUIDS IN TANK. DO NOT USE SOLVENTS IN THE CLEANER. DO NOT USE UNIT WITH GAS OR EXPLOSIVE CHEMICALS. DO NOT USE THE UNIT WITH DISINFECTANTS. KEEP THE AREA AROUND THE CLEANER CLEAN AND DRY.

In case of overflow, remove the water because of the electrocution risk. Clean up any spilt liquid. Water plus high voltages can cause electrical discharges.

NEVER PUT YOUR HANDS IN THE TANK WHEN THE UNIT IS RUNNING. NEVER DISASSEMBLE THE ULTRASONIC CLEANER. THE HIGH VOLTAGES INSIDE THE UNIT ARE HAZARDOUS.

If it is not possible to disconnect the power supply, disconnect the mains supply. If this is distant or not visible by the person carrying out the maintenance work, turn the switch to «OFF» and secure it with a padlock or, alternatively, lock the electrical cabinet and keep the key; always make sure the voltage is effectively absent.

DO NOT USE WITH A LIQUID LEVEL BELOW THE MINIMUM LEVEL OR RISE ABOVE THE MAXIMUM LEVEL INDICATED ON THE TANK. DO NOT USE WITH A LIQUID LEVEL BELOW THE MINIMUM LEVEL OR RISE ABOVE THE MAXIMUM LEVEL INDICATED ON THE TANK. DO NOT REMOVE THE LABEL AFFIXED TO THE UNIT. MAKE SURE THAT THE UNIT IS SUPPLIED WITH THE CORRECT VOLTAGE.

USE ORIGINAL SPARE PARTS ONLY.

If the equipment fails to operate properly, refer to the «Trouble Shooting» section of the manual for possible causes. If you should require further information or repairs, contact your supplier's specialised staff or the manufacturer technical office directly.

Warning

- The unit is built to the electrical safety standards set by official Standards Institutes. It is fitted with a two-pole plug to ensure complete earthing
- An essential safety requirement is to check that the electrical circuit has an efficient earth connection and that the circuit's current-carrying capacity and sockets match the unit's power rating, as stated on the specifications plate. Have this checked by qualified personnel.

FAILURE TO COMPLY WITH THE ABOVE POINTS ABSOLVES THE MANUFACTURER FROM ALL LIABILITY.

2.3 Electrical connection

Check that the supply voltage stated on the plate on the back of the unit corresponds to the voltage in the place of installation. The unit's power lead must be connected to a wiring circuit with an efficient earth connection, in compliance with binding regulations in the country of installation.

Caution: Do not bend the power lead or position it so that it is resting on objects of any kind.

CHAPTER 3

3.1 Description

The series of ultrasound tanks are the result of long and detailed studies into the use of this kind of washing system on an industrial scale. In particular, high-power ultrasound tanks have been developed for the dental market, for cleaning tools fouled by plaster, dentin powders, cement, etc. from tools and other objects. The stainless steel tanks are fitted with a high power ultrasonic transducer that converts the electrical power generated by an electronic tank circuit into mechanical ultrasonic energy, capable of cleaning the objects placed in the tank.



3.2 Technical specifications

Supply voltage	230 Volt
Number of phases	1
Supply frequency	50/60 Hz
Power input	80 W
Range of environmental conditions in which the unit was designed to operate	<ul style="list-style-type: none"> • Indoor use • Altitude up to 2000 m • Temperature: +5 ÷ +40°C • Max. relative humidity 80% for temperatures up to 31°C with a linear decrease to 50% at a temperature of 40°C • Max. variation in mains voltage ± 10% • Installation category (overvoltage category) II. Pollution level 1
Isolation class	I / IPX2
Fuses	2 - 5 x 20 - 2.5 A fuses
Weight	3,8 Kg
Transducer stack	1 with total power of 80 Watts
Operating frequency	30 - 32 KHz
Operating time control	0 - 30 min.
Heating element	---
Operating temperature	---
Colour	9016
Switch	Two-pole with protective cover
No aeration inlets, to prevent possible infiltration of fluids into unit.	

3.3 Ultrasound tank accessories

A range of supplementary accessories is available with the ultrasound tank. To ensure efficient operation, you are advised to procure at least the basket and cover together with the unit.

NAME	QUANTITY
Stainless steel basket, small mesh (180 x 100 x depth 60 mm)	1
Non-deformable transparent Plexiglas cover	1


3.4 Advantages of ultrasound cleaning compared to manual cleaning

- 1 Rapidly and thoroughly cleans surfaces which are intricate and difficult to reach: diamond-pointed cutters, endodontic files, manual Instrument, scissors, prostheses (crowns, inlays, removable prostheses), removable orthodontic appliances, laboratory instruments.
- 2 Prolongs the life of sharp instruments, thanks to the capacity of the ultrasonic apparatus to remove any residue from sharp surfaces.
- 3 Ultrasound cleaning is more rapid and efficient than manual cleaning, and is capable of removing any residue which interferes with chemical reactions.
- 4 Eliminates possible accidental cuts due to instrument brushing.
- 5 Eliminates the nebulization of dead organisms, generated by instrument brushing.


3.5 How ultrasound works

Ultrasound is a sound transmitted at a frequency out of the human hearing range (20 to 50 kHz); sound vibrations below 20 kHz are irritating and dangerous.


The radiation of waves causes high and low pressure stages in the liquid.

The ultrasound basin has a new vibration system called: SWEEP MODE, that is, vibration with change of frequency. With this type of operation, the stationary ultrasound waves are «broken up». This ensures greater uniformity of vibration and better and faster cleaning of the objects.

In the low pressure stage, millions of small bubbles are formed, that grow in size. This process is known as cavitation, i.e. the formation of cavities.

In the high pressure stage, the bubbles burst or implode, releasing enormous amounts of energy. These implosions react like a small army of brushes. They work in all directions, attacking all surfaces and entering all recesses and openings.

Ultrasound cleaning is very effective on hard materials and less so on soft and porous substances . The harder the surface, the more effective the ultrasonic «brushing» action. As metal, glass, ceramic materials and hard plastics are good conductors of sound, they are also the ideal materials for ultrasonic cleaning.

3.6 Ultrasonic cleaning methods

DIRECT METHOD

1. hanging
2. in basket

- Pour the cleaning solution into the tank
- The soiled articles are placed in the basket or hung

ADVANTAGES:


- Easy to perform
- Effective cleaning

INCONVÉNIENTS:

- The dirt removed remains in the tank
- Only one solution at a time can be used
- Highly acidic or caustic solutions can cause corrosion on the tank surface

With basket or hanging article

INDIRECT METHOD

1. beaker support
2. beakers

- Fill the tank with water plus a cleaning solution to improve cavitation
- The cleaning solutions are poured into the beakers or into the basket

ADVANTAGES:

- The dirt remains in the beaker
- Different solutions can be used at the same time
- E.g. One beaker with a cleaning solution, one beaker with a rinse solution

NEVER PLACE THE ARTICLES TO BE CLEANED ON THE BOTTOM OF THE TANK

DEGAS THE LIQUIDS EACH TIME THEY ARE CHANGED

3.7 Hints on using the unit

- Disconnect the unit from the electrical current before filling or emptying the tank.
- Fill the tank with lukewarm tap water, adding a cleaning solution to increase the cavitation.
- The operator is responsible for choosing the right cleaning agents for the articles or materials to be cleaned.
- If cleaning an article of value, you are advised to do a trial run with a small test sample and then proceed with the remaining material.
- Always use baskets or beakers; if these are not used, the articles must be hung.
- Never let articles rest on the tank bottom.
- Always keep the level of the liquid between the MIN and MAX levels indicated on the tank.
- Clean the tank at regular intervals and change the cleaning solution frequently to maintain the effectiveness of the cleaning action.
- Do not restrict the ventilation, which must be kept at a suitable level.
- Do not use the tank in very dusty surroundings. Accumulated dust could impair the efficient operation of the unit.
- Degas every time the liquids are changed, running the unit with the new solution for about 5 /endash 10 minutes before putting any articles in the tank.
- Rinse the cleaned article in drinkable water when it has completed its cleaning cycle (demineralised water is preferable, to avoid staining), then dry it thoroughly. If the article is a medical instrument, put it through the sterilisation cycle.

OPERATING THE APPLIANCE WITH THE ARTICLES PLACED ON THE BOTTOM OR WITH THE LIQUID LEVEL OUT OF THE RANGE INDICATED ON THE TANK WITH THE MIN AND MAX MARKS CAN IRREVERSIBLY COMPROMISE THE CLEANER.

CHAPTER 4

4.1 Instructions for use

THE APPLIANCE MUST NEVER BE ALLOWED TO WORK WITHOUT WATER.

THE APPLIANCE MUST NEVER BE ALLOWED TO WORK WITHOUT THE COVER ON.

- 1 Fill the tank with lukewarm tap water (max 60°).
- 2 Add the cleaning solution to the water. N.B.: the level of the liquid must lie between the MIN and MAX marks indicated on the tank, when all the articles with the relative baskets/cups, have been immersed.
- 3 Switch the unit on with the luminous ON/OFF switch on the front, which lights up when the power is on.
- 4 Set the operating time with the timer to start generating ultrasounds. To do this, turn the central knob until the notch points to the required time. The blue light then shines as well as the green one. The generating of ultrasounds stops automatically when the timer has elapsed and returns to zero. To switch the unit off completely, turn the main switch on the back to OFF.

Note:

- Degas every time the liquids are changed, running the unit with the new solution for about 5 /endash 10 minutes before putting any articles in the tank.
- For cleaning blood, you are advised to keep the cleaning solution at a low temperature (30-40°C) and then switch the unit off after use.
- For cleaning cements, alginates and calcareous substances, take the solution to a high temperature.

4.2 Definition of the colours front panel

FRONT PANEL COLOURS		MEANING
Green light shining		Rear switch ON (power ON)
Green light and blue light shining		Ultrasound working
Red light shining		Overheating safety trip. Turn the unit off and on again using the rear switch.

CHAPTER 5

5.1 Routine maintenance and precautions

WEAR SUITABLE PPE DURING MAINTENANCE.

- Each time the liquids are changed they must be DEGASSED. If this is not done, at least for the first cleaning cycle, the cleaning power of the ultrasound tank will be reduced.
- Change the cleaning solutions at regular intervals. Empty the tank and dry the bottom with a cloth. Do not allow the impurities removed from the cleaned instruments to be deposited on the bottom of the tank in a significant quantity for an extended period of time. This will reduce the ultrasound cleaning power. Also, after a certain time, the liquids become depleted and no longer have any cleaning effect. It cannot be known beforehand how long the liquids will last, as this depends on the number and duration of the cleaning cycles.
- Always comply with binding regulations in the country of use when disposing of depleted liquids.

Note: You are advised not to use liquids with a high acid content, or at a temperature of over 80°C. In fact, the closer the liquid to boiling point, the less effective the ultrasound cleaning (as water comes to the boil, it absorbs air, which impedes the cleaning action of the ultrasound wave.).

- You are advised to use the cover during the work cycle, to contain the vapours produced.
- Wipe the exterior of the unit with a slightly dampened cloth. Before doing this, make sure that the tank is disconnected from the electrical power supply.
- All maintenance, including the changing of the power supply lead, must be performed by authorised, qualified personnel.

5.2 Problem solving

PROBLEM	CAUSE	REMEDY
1. Unit does not work and luminous switch does not light up	Power not on	Mains power failure
	Switch broken	Change switch
2. Unit does not work and luminous switch remains lit	Safety fuses	Check for short circuits or faults, then replace the fuse with another one of the same grade.
3. Ultrasound generators do not work, with luminous switch lit up and heating working	Safety fuse on generator board	Check for short circuits or faults, then replace the fuse with another one of the same grade
	Faulty generator board	Change generator board
4. Ultrasound effectiveness reduced	Solution not degassed	Degas the liquids, following the instructions in the section «Hints on using the unit»
	Solution depleted	Change the solution, cleaning the tank with a cloth
	Solution not suitable	For cements, plasters, alloys etc., use the specified liquids to get the best results.
	Transducers working at wrong frequency	Contact Assistance Service for repairs

SWITCH OFF THE POWER BEFORE CARRYING OUT ANY WORK ON THE UNIT. HAVE THE WORK DONE BY SPECIALIST ENGINEERS ONLY

HIGH INTERNAL VOLTAGES

CHAPTER 6

6.1 Scrapping instructions

The ultrasound unit has been manufactured using ferrous materials, electronic components and plastic materials. No particular precautions are required when demolishing the unit.

Do not dump the unit in an unguarded place; consign it to a waste disposal agency for scrapping. Always refer to the binding laws in the country of use when scrapping equipment.

The  symbol on the appliance means that it must be disposed of as “sorted waste”.

The user must therefore send (or instruct other people to send) the unit to one of the sorted waste collection centres set up by the local councils, or send it to the dealer against the purchase of an equivalent unit (European Union only)

Sorted waste collection and the subsequent treatment, recovery and disposal operations facilitate the production of equipment using recycled materials and limit the negative environmental and health effects that may be caused by improper waste management.

Abusive dumping by the user will be punished according to law.

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Distributed by CADENCE
2 bis, chemin du Loup
93290 Tremblay-en-France - France

WARRANTY STATEMENT

This instrument described below has been fully inspected and conforms to the current products specification.
This instrument is guaranteed for its designated use, against original defects in materials and workmanship for a period of 24 months from date of purchase.

Repairs due to misuse or disassembling of the equipment, or by the parts not provided by the manufacturer, are not covered by this warranty.

The following parts are not included in warranty: Stainless steel basket, plexiglas cover.

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Distributed by CADENCE
2 bis, chemin du Loup
93290 Tremblay-en-France - France

(tear off here)

INFORMATION CARD

(to be returned)

Instrument Type:

Serial N°:

Inspector:

Power: Power:

Customer name:

Adresse:

Purchase date:

Please fill this « Information card » and send it to the sales representative for user's record and future service.

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