



Choice of Millions of Dentists Worldwide



Guilin Woodpecker Medical Instrument Co., Ltd.
Information Industrial Park, National High-tech
Zone, Guilin, Guangxi, 541004 P. R. China

Registrant/manufacturer/after-sales service provider:
Guilin Woodpecker Medical Instrument Co., Ltd.
Tel:
Sales Dept.: 0773-5873196/ 2350599
Fax: 0773-5822450
After-sales Service Dept: 0773-5827898/13978361362
E-mail: woodpecker4@glwoodpecker.com
Website <http://www.glwoodpecker.com>



Upgraded star product

Brushless Endo Motor with
two-in-one function and compact
handpiece

Endo Radar Pro

- Two-in-one function enables the root canal preparation with simultaneous length control
- Super-mini handpiece
- Brushless motor
- Cordless



ENDO RADAR
PRO


SYSTEM
W3-Pro WOODPECKER
TORQUE 1.5 Nm
SPEED 350 RPM

Auto Rev
Rev
EnAp
SET
SYSTEM
FILE
TORQUE
SPEED
ANGLE


WOODPECKER
Endo Radar Pro

Two-in-one function

- ① Root canal preparation with simultaneous length control



- ② Length determination only



A. Auto start

Automatically start when the file is inserted into the canal and stop when the file is taken out.

B. Auto apical reverse and auto apical stop

Automatically reverse and stop when the file reaches the apical stop.

C. Auto apical slowdown

The file slows down as it gradually approaches the apical stop.





Endo Radar

Ultra-compact and cordless handpiece enables better hand feeling


Handpiece diameter reduced by 2mm, circumference reduced by 7mm
Control the handpiece with wireless signal, without the restraint of tail cord

Endo Radar Pro


Diameter reduced by 2mm
Circumference reduced by 7mm

Super-mini contra-angle

The contra-angle is only 8mm in diameter and 9.7mm in height



Head and neck of ordinary contra-angle



Head and neck of Woodpecker super-mini contra-angle



6:1 speed ratio


Ensure low speed and stable rotation under high torque.



360° rotatable handpiece

The contra-angle can be rotated directly without being pulled out, which is convenient for the preparation of different teeth.






High performance brushless motor

Adjustable reciprocating angle with 10° interval
 Adjustment range of forward and reverse reciprocating angle: 20° - 340°
ATR mode: After the preset torque is reached, the file will reciprocate according to the set angle.
 Accurate control of the rotation angle reduces the risk of instrument separation and enables more efficient and safer preparation.

Smoother rotation, low noise and low vibration



High-performance brushless motor with high torque, high speed and fast response
 Speed: 100-1200rpm Torque: 0.4-5.0 N.cm

High-performance brushless motor for low wear, low heat and longer service life

Battery life:

Battery life of brushless motor = battery life of brush motor x10





All touch button design HD LCD screen

- Touchable smart panel enables comfortable hand feeling and easy cleaning
- The luminous panel guarantees easy identification of buttons


Wireless charging for longer endurance

- Built-in 1200mAh lithium battery in the handpiece enables wireless charging
- 2000mAh large capacity battery built in the base
- The handpiece can be charged when placed on the main unit, enabling easy use and long durability.
- The use of the machine is not affected by power failure.
- The brightest Endo Motor in the dark

Built-in 11 user-defined modes to adapt to mainstream file systems

Easy to choose and easy to operate

Operation steps (take W3-Pro file as an example)



1. Use # 10 stainless steel K file to probe and negotiate the root canal, and conduct preliminary determination of working length;
2. Use cavity access file 17/.12 to pre-enlarge the coronal 1/3 of the root canal and remove the ferrule at the root canal orifice;
3. Determine the working length and confirm the patency of the root canal;
4. Use glide path file 12/. vt to negotiate and enlarge root canal access;
5. Use glide path file 16/.vt to negotiate and enlarge root canal access;
6. Use file 18/.05 to negotiate the root canal to the working length by means of "brushing" and "pecking";
7. Use file 25/ .06 to shape the root canal to the working length by means of "brushing" and "pecking";
8. Use file 35/.04 to shape the root canal to the working length by means of "brushing" and "pecking".

Consciously use the "brushing" method to enlarge the lateral space of the root canal to ensure safer root apex preparation; consciously use the "pecking" method to avoid apical deviation of root canals more effectively.

