BiodentineTM XP

Active Biosilicate TechnologyTM

fr	Substitut dentinaire bioactif	1
en	Bioactive dentine substitute	4
es	Sustituto bioactivo de la dentina	7
de	Bioaktiver Dentinersatz	10
it	Sostituto bioattivo della dentina	13
pt	Substituto bioativo de dentina	17
nl	Bioactief dentinesubstituut	20
pl	Bioaktywny substytut zębiny	23
el	Βιοενεργό υποκατάστατο οδοντίνης	26
bg	БЮАКТИВЕН ЗАМЕСТИТЕЛ НА ДЕНТИН	29
cs	Bioaktivní náhrada dentinu	33
da	Bioaktiv dentinerstatning	36
et	Bioaktiivne dentiini aseaine	39
hu	Bioaktív dentin helyettesítő	42
lt	Biologiškai aktyvus dentino pakaitalas	45
lv	Bioloģiski aktīvs dentīna aizstājējs	48
ro	Substitut de dentină bioactiv	51
sk	Bioaktívna náhrada dentínu	54
sl	Bioaktivni nadomestek dentina	57
sv	Bioaktiv dentinersättning	60
hr	Bioaktivna zamjena za dentin	63
no	Bioaktiv dentinerstatning	66

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CH REP Mandataire CH

en

Biodentine™ XP is a bioactive dentine substitute based on the "Active Biosilicate Technology™" innovation.

PRESENTATION

Biodentine™ XP 500 (red piston): box of 10 foils, each foil contains 1 cartridge of 700 mg of powder and 190 µL of Liquid.

Biodentine™ XP 200 (black piston): box of 10 foils, each foil contains 1 cartridge of 468 mg of powder and 133 µL of Liquid.

COMPOSITION

The cartridge contains both the powder and liquid components:

- The powder component is composed of tricalcium silicate, zirconium oxide, calcium oxide, calcium carbonate and iron oxides.
- The liquid component is composed of calcium chloride, polycarboxylate and purified water.

INTENDED USE

Bioactive dentine substitute.

TARGET GROUP(S)

Target population

Biodentine™ XP is intended to be applied on temporary teeth (from 2 years old), on permanent teeth (mature or immature) or on non-vital immature permanent teeth depending on indication. Please refer to the dedicated section (§Indications for use).

Particular population:

Limited clinical data is available for patient population aged from 2 to 6 years old.

No clinical data is available for pregnant or breastfeeding women.

Intended users

Dental healthcare professionals only.

CLINICAL BENEFIT

- Preservation of tooth vitality (except for apexification, perforations and furcations)
- Preservation of tooth functionality
- Revitalization of necrotic teeth

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A Summary of Safety and Clinical Performance is available, using Basic UDI-DI Code 37601407000004G9, in the European database on medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed> or upon request (sscp@septodont.com) until EUDAMED is fully functional.

PERFORMANCE CHARACTERISTICS

Technical performances / specifications

- Compressive strength:

At 24h, compressive strength \geq 150MPa

After 28days, compressive strength: between 260MPa and 300MPa which is similar to the sound dentine.

- Sealing ability: creation of a tight seal within dentinal tubules.
- Setting time: the average setting time of the product is 12 minutes, starting from the end of mixing. The setting time may vary (from 9 to 25 minutes) depending on some conditions, including the storage time and conditions.
- Working time: \geq 1 min

Biological performance

The bioactivity of Biodentine™ XP relies on hydroxyapatite formation promoting a favourable environment for angiogenesis and osteogenesis, isolating the root canal from the surrounding tissues and stimulating the healing process of damaged apical tissues.

Clinical performance

- Preservation of the pulp
- Preservation of the amelodentinal structure
- Preservation of the root canal structure
- Root edification
- Sealing ability necessary to create a favourable environment required for the Revitalization procedure

INDICATIONS FOR USE

Biodentine™ XP 500 and Biodentine™ XP 200 are used for the same indications, the only difference is quantity.

In the crown:

For temporary teeth (children from 2 years of age) and permanent teeth (immature or mature):

- Permanent dentine restoration under composites or Inlay/Onlay
- Temporary dentine-enamel restoration
- Restoration of deep and/or large coronal carious lesions (sandwich technique)
- Restoration of cervical or radicular lesions

On the pulp:

For temporary teeth (children from 2 years of age) and permanent teeth (immature or mature):

- Pulp capping (direct and indirect)
- Pulpotomy for diagnosed symptoms of reversible pulpitis and irreversible pulpitis where bleeding is controlled within 5 minutes.

In the root:

For permanent teeth (immature or mature):

- Repair of root perforations
- Repair of furcation perforations
- Repair of perforating internal resorptions
- Repair of external resorption
- Root-end filling in endodontic surgery (retrograde filling)

For immature permanent teeth:

- Apexification
- Revitalization procedure by means of revascularization of permanent immature teeth with necrotic pulp

CONTRA-INDICATIONS

- Hypersensitivity to one of the ingredients.
- Restoration of large loss of tooth substance subjected to high stresses.
- Pulp space that requires post or core for the final restoration.

EXPECTED LIFETIME

At least 10 years.

INSTRUCTIONS FOR USE

Refer to Precautions for use 1 and 2

 Scan the QR code on the Protocol Card (Document "PROTOCOL CARD") to watch the video.

1. Activation of the cartridge with the Biodentine™ Mixer

- 1.1. Open the lid of the Biodentine™ Mixer.
- 1.2. Put and align the Biodentine™ XP cartridge in the activation socket (see Protocol Card step 3).
- 1.3. Lock the cartridge manually by making a clockwise quarter turn.
- 1.4. Manually pull the lever down all the way (until it stops), you may hear a slight click (see Protocol Card step 4).
- 1.5. Put the lever back in its initial position.
- 1.6. Unlock the cartridge from the activation socket by making an anticlockwise quarter turn (See Protocol Card step 7).
- 1.7. Check if the cartridge is properly activated: the cartridge piston should be pushed inside. Refer to Precautions for use 4.
- 1.8. Once the cartridge is activated, proceed immediately to Mixing. Refer to Warning 9.

2. Mixing

- 2.1. Transfer the cartridge from the activation socket to the Mixing arm of the Biodentine™ Mixer (see Protocol Card steps 7 and 8).
- 2.2. Lock the cartridge manually by making a **clockwise quarter turn**.
- 2.3. Close the lid of the Biodentine™ Mixer.
- 2.4. Press the START button on the Biodentine™ Mixer, (see Protocol Card step 9).
- 2.5. The START button will progressively become green indicating the elapsed mixing time (mixing time: **30 seconds**).

- 2.6. When mixing is complete, the START button will blink green and beep twice.
- 2.7. Open the lid of the Biodentine™ Mixer.
- 2.8. Unlock the cartridge from the mixing arm by making an **anticlockwise** quarter turn, then remove it.
- 2.9. Close the lid of the Biodentine™ Mixer. Turn off the mixer when further operation is no longer needed.

3. Insertion of a Biodentine™ cartridge into the Biodentine™ Gun

- 3.1. Immediately after mixing, raise the Biodentine™ Gun lever until it stops in a horizontal position.
- 3.2. Align the Biodentine™ XP Cartridge in the Biodentine™ Gun (see Protocol Card step 10).
- 3.3. Insert the cartridge into the Gun and lock the cartridge making a **clockwise quarter turn** (see Protocol Card step 11).

4. Priming

- 4.1. Hold the Biodentine™ Gun, the lever in one hand and the handle in the other, (see Protocol Card step 12) and **progressively** push the lever down. **Refer to Warning 11.**
- 4.2. Stop priming when the dentine substitute appears in the nozzle (see Protocol Card step 13).

5. Extrusion

- 5.1. Hold the Biodentine™ Gun with one hand (see Protocol Card step 14).
- 5.2. With some effort, the cartridge nozzle can be bent manually (due to the thickness of the nozzle) to facilitate the placement of the dentine substitute (see Protocol Card step 15) **Refer to Warning 10.**
- 5.3. Progressively push down the lever to allow the dentine substitute to flow out of the cartridge. A homogeneous and creamy dentine substitute is obtained. **Refer to Warning 4.**

6. Application

Apply the dentine substitute directly into the affected tooth cavity with the Biodentine™ Gun. (see Protocol Card step 16). Alternatively, you may apply Biodentine™ XP with an amalgam carrier, a spatula, or a Root Canal Messing Gun.

7. Cartridge removal

Release the lever and unlock the cartridge from the Biodentine™ Gun by making an **anticlockwise** quarter turn and then removing the cartridge (see Protocol Card step 17).

8. Disposal

Dispose the cartridge in an appropriate bin. Refer to the Disposal section.

WARNINGS AND PRECAUTIONS FOR USE

Precautions for use

1. Wear the appropriate personal protective equipment while using the device because it may cause eye irritation.
2. All the instructions for use are intended for application within the working area, away from the patient.
3. Water contamination slows down the setting of the material. Prevent exposure to water and fluids during the initial setting stage.
4. If the cartridge is not correctly activated, it may fall from the mixing arm.
5. Do not prime a cartridge that has not been locked properly with a **clockwise quarter turn** in the Biodentine™ Gun.
6. Refer to the Biodentine™ Mixer and Biodentine™ Gun instructions respectively for use for the Biodentine™ Mixer and Biodentine™ Gun preparation.

Warnings

1. Use a rubber dam to ensure that the dentine substitute does not come into contact with the mucosa and is not swallowed.
2. Single-use product: to avoid cross contamination, do not use the mixed dentine substitute on several patients.
3. Biodentine™ XP is only intended to be mixed with the Biodentine™ Mixer and extruded with the Biodentine™ Gun. Otherwise, its technical specifications are not guaranteed.
4. If the consistency of the obtained dentine substitute is other than homogeneous and creamy, do not use it.
5. Wait until the end of the setting time before performing permanent enamel restoration or removing the matrix (when used).
6. For revitalization procedure, wait until the formation of a blood clot before applying Biodentine™ XP. Anaesthetics with vasoconstrictors may prevent blood clot formation.
7. For endodontic procedures, flush with a sterile saline solution to remove any irrigant residue and dry without totally dehydrating the root canal.
8. If the cartridge falls on the floor, clean it with ethanol to avoid contamination.
9. If **3 minutes pass** between cartridge activation and mixing, do not use the cartridge, otherwise the technical performance of the dentine substitute is not guaranteed.
10. Do not bend the nozzle to 90° more than 4 times to avoid breaking the nozzle.
11. During the priming step, do not push the lever down too suddenly; push it progressively to avoid unintended extrusion of the dentine substitute.

SIDE-EFFECTS

Possible hypersensitivity to Biodentine™ XP or ingredients (*estimated frequency: less than 1/100 000*).

Please report any serious incident that occurs involving the device to Septodont (vigilance@septodont.com) and to the local competent authority for medical devices. Please inform the patient of potential side-effects related to Biodentine™ XP.

STORAGE

Store between 5°C and 25°C and protect from moisture.

Any remaining material shall be discarded.

DISPOSAL

Dispose of Biodentine™ XP in an appropriate bin if it is contaminated, otherwise it can be disposed of in a regular bin.

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