# SS\_ONE SS\_KISS

MANUALE D'USO\_ITA  
USER MANUAL\_ENG  
MANUAL DE INSTRUCCIONES\_ESP

MANUALE D'USO	p. 1
USER MANUAL – ENG	p.33
MANUAL DE INSTRUCCIONES – ESP	p.65
CERTIFICATI _ CERTIFICATES _ CERTIFICADOS	p.97
GARANZIA _ WARRANTY _ GARANTÍA	p.104
CERTIFICAZIONE DI CONTROLLO CONTROL CERTIFICATION _ CERTIFICACIÓN DE CONTROL	p.106
REGISTRO VERIFICHE ELETTRICHE REGISTER OF ELECTRICAL TESTS REGISTRO DE PRUEBAS ELÉCTRICAS	p.108
ASSISTENZA _ ASSISTANCE _ ASISTENCIA	p.110
SCHEDA DI INSTALLAZIONE _ INSTALLATION SHEET FICHA DE INSTALACIÓN	p.112

# ITALIANO

MANUALE D'USO

# Indice

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## USER MANUAL

# Index

<b>1 _ Safety information</b>	<b>36</b>
<b>2 _ Technical Characteristics of the Product</b>	<b>38</b>
<b>3 _ Product Introduction</b>	<b>39</b>
3.1 _ Product Introduction	40
<b>4 _ Basic Operations</b>	<b>44</b>
4.1 _ The Dental Unit	44
4.1.1 _ Starting	44
4.1.2 _ Configuration of the Fuse into the Floor Unit	45
4.1.3 _ Footswitch	46
4.1.4 _ The Headrest of the Unit	47
4.2 _ The Operator Unit	48
4.2.1 _ Handpieces Activation	48
4.2.2 _ Flexible Arm Brake	48
4.2.3 _ Main Control Panel	49
4.2.4 _ Foot Control	49
4.3 _ Water Unit	50
4.3.1 _ Water Supply System	50
4.4 _ Assistant Unit	51
4.5 _ The Cuspidor	52
<b>5 _ Control Panel and Settings</b>	<b>52</b>
5.1 _ Control Panel	52
5.1.1 _ Chair Positions	54
5.1.2 _ Boiler for hot water into the Water Unit	55
5.2 _ Settings	55
5.2.1 _ Chair Setting	56
5.2.2 _ Cuspidor Setting	56

<b>6 _ Adjustment</b>	<b>57</b>
6.1 _ Doctor's Unit	57
6.1.1 _ Handpieces Adjustment	57
6.2 _ The Cuspidor	59
6.2.1 _ Adjustment of the Cuspidor Rinse and Filling 59 of the Glass	
<b>7 _ Cleaning and Disinfection</b>	<b>60</b>
7.1 _ Cleaning and Care on the Cordons of the Handpieces	60
7.2 _ Cleaning of the Suction Cannulas	61
7.3 _ Oil Collector	61
7.4 _ Cleaning of the Suction Filter	62
7.5 _ Cleaning of the Cuspidor	62
7.6 _ Cleaning of the Lamp and the Unit Surfaces	63

## 1 \_ Safety Information

Before using the unit, please read the manual carefully.

- Make sure that the unit is positioned firmly on the ground.
- This unit does not include those components which are not covered in the use. In case of problem, do not attempt to dismantle the unit by yourself.
- All electrical components and connections must be certified by a qualified technician.
- Only accessories admitted by Simple&Smart can be used.
- After each use of the unit, the patient chair must be moved to the lowest position and the water, the air and the power supply must be turned off.
- Do not use this unit in environments with high temperature, high humidity, in unstable conditions or in a dusty place.
- Only the authorized spare parts can be used.
- The unit must not come into contact with any object while the chair is moving.
- Do not touch the bulb and the surface of the reflector of the lamp while it's working.
- The device should not be used in environments where are located equipments that generate magnetic or electromagnetic fields.
- The equipment can not be used for purposes other than those intended.
- **The dental unit you purchased is a medical device subject to post-sales monitoring by Simple & Smart. In particular, it is required to promptly provide information regarding malfunctions that could compromise, or have already compromised, the safety and performance of the dental unit.**

### **Important Note on the manual and its translation.**

In case of disputes, it is considered as a reference manual, the Italian version, or the English one. Any other translation, made by third persons, will not be taken into consideration.

**These instructions describe how to properly use the dental unit SS\_ONE/  
SS\_KISS.**

The information, technical characteristics and the pictures, contained in this publication are not binding. SIMPLE&SMART srl reserves the right to make technical changes and improvements without modifying these instructions.

The manufacturer pursues a policy of constant improvement of its products, so, it is possible that some instructions, specifications and pictures in this manual may differ slightly from the product you purchased.

The manufacturer also reserves the right to make any change to this manual without notice. This equipment must be used only by personnel (medical and paramedical) adequately trained.

**SIMPLE&SMART s.r.l. decline any responsibility for any possible problem caused by disregard of the warnings and instructions contained in this manual.**

## 2 \_ Technical Characteristics of the Product

Table 1 \_ Technical characteristics

Power dissipation	AC220-230V,50/60Hz,350VA o AC110V~120V 50/60Hz, 350VA
Illumination power of operating lamp	8000-25000 Lux
Intensity of the X-ray viewer	$\geq 2000\text{cd/m}^2$
Carrying capacity of the unit	1323N (135kg)
Carrying capacity of the main tray arm	2.5kg
Maximum air power to the turbine	Idle speed $\geq 35 \times 10^3$ r/min (Air pressure: 0.22Mpa)
Minimum air power to the handpieces	Idle speed $\geq 18 \times 10^3$ r/min (Air pressure: 0.3Mpa)
Highest position of the sitting patient from the ground	800mm
Lowest position of the sitting patient from the ground	400mm
Backrest movement range	$109^\circ - 172^\circ$
Headrest extension range	0-100mm
Air supply pressure	0.55-0.80Mpa
Water supply pressure	0.2-0.4Mpa
Main line fuse	250V 6A( 220V input) 250V 10A( 110V input)
Power cable	227 IEC(RVV)3x0.75mm <sup>2</sup>

## 3 \_ Product Introduction

Thank you for choosing a Simple&Smart equipment.

This guide contains basic operating instructions, instructions for the setting of the control panel, adjustments and instructions about maintenance of SS\_ONE / SS\_KISS systems. This includes:

- The chair
- Assistant unit
- Water unit, including the cuspidor
- The operating lamp
- The doctor unit with hanging type or top mounted cordons
- Foot control

Any optional or different configurations can be found onto the catalog or online at [www.simplesmart.it](http://www.simplesmart.it)

### 3.1 \_ Product Introduction

Figure. 1 \_ SS\_ONE/SS\_KISS Unit


Table 2 \_ SS\_ONE/SS\_KISS Configuration

Position	Description
A	Chair with movement device
B	Assistant unit
C	Cuspidor
D	Lamp (Halogen or LED)
E	Water unit
F	Main unit arm with Doctor unit (hanging type or top mounted cordons)
G	Endoscope system (Monitor e camera)
H	Floor Unit with integrated controls
I	Multi Function Foot control (or FARO pedal)

Figure 2 – Main Unit, Hanging Type tray

Table 3 \_ Hanging type configuration

Position	FUNCTION
A	Doctor Tray
B	Main unit control panel
C	Handle with brake
D	Support of the arm
E	Instruments support
F	Instruments

Figure 3 – Main Unit, Top Mounted tray

Table 4 \_Top mounted configuration

Position	FUNCTION
A	Doctor Tray
B	Main unit control panel
C	Handle with brake
D	Support of the arm
E	Support of the instruments
F	Instruments

Figure 4 – Assistant Unit


Table 5 \_ Assistant tray configuration

Position	FUNCTION
A	Strong suction cannula
B	Weak suction cannula
C	3 ways Syringe
D	Support of the tray
E	Intraoral camera
F	Empty space for any optional (curing light)

## 4 \_ Basic Operations

This section includes basic operating instructions for the use of the chair, the main unit, the assistant unit and the water unit with the cuspidor.



For the instructions concerning the control panel, see the relevant chapter.

### 4.1 \_The Dental Unit

#### 4.1.1 \_ Starting

Connect the dental unit to the electrical socket using the special plug and turn on the device through the power button located inside the base. The activation of water and air is done by means of the lever placed outside the base. It is strongly recommended to reposition the lever in the OFF position at the end of the operation.


Figure 5 \_ Floor Unit with Integrated controls

Activation / Switch off lever for hydraulic and pneumatic system.

## 4.1.2 \_ Configuration of the Fuse into the Floor Unit

Figure 6 – Configuration of the fuse into the floor unit


- F1 e F2 \_ (220V AC input) Fuse 6A
- F3 \_ 24V AC (yellow) Fuse 10A
- F4 \_ 24V AC (white) Fuse 10A
- F5 \_ 0V (black) Fuse 10°
- F6 \_ empty

## 4.1.3 \_ Footswitch

The movement of the chair can be done using the pedal switch, which allows to move the seat and the backrest of the chair up and down.

Figure 7 \_ The pedal switch

The backrest movement consists in a variation of angle toward the floor, ranging from 8 to 71 degrees. In some configurations the switch is not present, in some others, it is positioned on the rear plastic cover instead of on the pedal.


The unit is not equipped with an anti-crushing safety system, therefore, avoid to leave any component beneath the seat and pay close attention to those parts of the body that could be crushed during the movement of the seat.

## 4.1.4 \_ The Headrest of the Unit

Figure 8 – Adjustment of the headrest

- To adjust the height and the angle of the headrest, loosen the block counterclockwise.
- Push and press the headrest to reach the desired height and the angle desired.
- Close the block by tightening clockwise.

- Do not draw the headrest for more than 100 mm otherwise you will remove it from its housing.
- After adjusting the headrest, it shall be locked firmly to avoid accidents.

## 4.2 \_ The Operator Unit


### 4.2.1 \_ Handpieces activation

The handpieces automatically activate when you pull forward the cords and you lift them from their rest position. By pressing the foot control you can operate with the handpiece previously activated.

### 4.2.2 \_ Flexible Arm Brake

The flexible arm brake allows you to provide various weight settings (2.5kg. maximum) on the main unit tray, keeping the position. The brake limits the vertical movement but allows the horizontal positioning of the control head.

Figure 9 \_ Switch of the brake on the unit

## 4.2.3 \_ Main Control Panel

The control panel allows you to control the chair, the cuspidor, the lamp and the X-ray film viewer through the relative buttons.

For complete instructions on the control panel, see the section on the following pages.

## 4.2.4 \_ Foot Control

The foot control modulates the incoming air to the active handpiece and provides a signal which activates the air and water flow. The control is operated by a light foot pressure applied to any part of the pedal disc. The foot control is also equipped with a switch for the water activation and a button for the chip air. In case the unit is equipped with a FARO pedal, it is activated through an horizontal movement of the lever and the water activation through a vertical movement. The FARO pedal is not equipped with further controls.

- To turn off the water supply to the tools, you must move the lever of the relative switch.
- To turn on the water supply, you have to move the switch lever in the opposite direction.

Figure 10 \_ Multi function pedal

## 4.3 \_ Water Unit

### 4.3.1 \_ Water Supply System


The 2 independent tanks, one of 600 ml and the other of 2 L capacity, can supply water to the handpieces and syringes. Fill the tank mounted outside the water unit with the desired solution and activate it using the appropriate switch located next to the tank itself (In case you want to activate the tank for water solutions) or fill the tank inside the water unit and activate it by means of the special switch located near the tank (in case you want to use the disinfection tank). The use of the two tanks replaces the use of the standard water supply network.

Figure 11 \_ Water supply system

#### Tank 1 \_ for water based solutions

Switch for activation of tank 1

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Switch for activation of tank 2 (inside the water unit)

- The default setting (both switches in OFF position) ensures the supply from the urban water network, by activating the respective switches you will start the water supply from the tank.
- At the end of the operation, do not forget to turn off the switches to avoid entering of undesired solution inside the patient's mouth.

## 4.4 \_ Assistant Unit

The assistant unit includes: the control panel, the strong suction cannula, the weak suction cannula and the 3-ways syringe. The suction cannulas are connected to the water unit through 2 airtight filters, for the separation of the solid parts from the liquid ones, one for each cannula.

Figure 12 \_ Assistant unit

It is advisable to clean the filters every day. We also remind you to take note of the procedures for disinfection and sanitization of the chair listed in the relevant section.

## 4.5 \_ The Cuspidor

The cuspidor can be activated and programmed from the control panel.

The cuspidor is activated by pressing the button 

while the glass filling is activated by pressing 

The glass filling is preset by the production with an automatic timing of 3 seconds and the spitton rinsing with a timing of 12 seconds. Both the glass filling and the spitton rinsing can be setted again according to your needs. For the operation and the programming, see the relative section on the following pages.

## 5 \_ Control panel and settings

### 5.1 \_ Control panel

The control panel is operative when the green LED light is turned on.


Figure 13 \_ Control panel of the top mounted main unit

The control panel can give commands for the movement of the chair, the control of the cuspidor, the operating light, the X-ray film viewer and the boiler for hot water. Moreover, it has the memorize some positions of the chair according to your needs.

For the programming, refer to the following pages.


Figure 14 \_ Control panel on the hanging type main unit

The control panel can give commands for the movement of the chair, the control of the cuspidor, the operating light, the X-ray film viewer and the boiler for hot water. Moreover, it has the memorize some positions of the chair according to your needs.

For the programming, refer to the following pages.

## Control panel of the assistant unit


Figure 15 \_ Control panel of the assistant unit. For both the models, hanging type or top mounted.

The control panel of the assistant unit, can give commands for the movement of the chair, the control of the cuspidor, the operating light and the boiler for hot water.

The programming is not possible from the assistant unit.

The shape of the control panel can vary depending on the version of the model.

## 5.1.1 \_ Chair Positions

Movement icons on the control panel allow you to manually move the seat of the chair and the backrest up and down.

Table 6 \_ Movement of the chair

Button	Function
	Backrest Up
	Backrest Down
	Seat Down
	Seat Up

## Buttons for the chair position programming

Table 7 – Control panel of the chair positions

Button	Function
	Automatically returns the chair to the position of entry and exit of the patient
	Automatically brings the chair to the rinse position. By pushing it again, it takes back the chair in the last working position.
	Program 1 _ chair position
	Program 2 _ chair position
	Program 3 _ chair position

## 5.1.2 \_ Boiler for Hot Water into the Water Unit

Press and hold the button  for a few seconds to activate the water heating function. The LED signal on the control panels indicates that the boiler is working. Hot water can be supplied for a few minutes, the Led remains lit up until the water cools down. You can press the button again to turn on the boiler or to switch of the heating function.


The green turned on LED, placed beside the icon of the boiler on the control panels, indicates that the boiler is working. It is recommended to turn off the water boiler when it is not in use.

## 5.2 \_ Settings

To activate the programming mode, press and hold the SET button until you turn on the corresponding LED indicator.

Figure 16 \_ LED light indicating the activation of the programming mode on the control panel

## 5.2.1 \_ Chair Setting


### Programming procedure

- Move the seat and the backrest to the working position you want to save.
- Press the "SET " key until the LED indicator on the control panel glows.
- Press the P1 or P2 or P3 key and then press the " SET " button again to finish the setting procedure. The LED indicator turns off and the P1 or P2 or P3 key is programmed with the selected position.
- Repeat the same operation to set the other positions.


## 5.2.2 \_ Cuspidor Setting

### Glass filling

- Press the "SET " key until the LED indicator on the control panel glows.
- Press the button for the water supply and keep it pressed until the water level reaches the desired amount, then release the button; the water supply will stop.
- Press the " SET " button again. The LED indicator turns off and the setting procedure is finished.

## Basin rinsing


- Press the "SET " key until the LED indicator glows.
- Press and hold the button for the basin rinse until you want.
- Press again the "SET" button to finish the setting procedure.

## 6 \_ Adjustment

### 6.1 \_ Doctor's Unit

#### 6.1.1 \_ handpieces adjustment

Figure 17 – Handpieces adjustment on the hanging type main unit


A : adjustment of the water/air spraying of all cords

B : regulation of the water flow of each single cord

C : regulation of the air flow of each single cord

Figure 18 – Handpieces adjustment on the top mounted main unit

A : adjustment of the water/air spraying of all cords

B : adjustment of the water flow of each single cord

C : adjustment of the air flow of each single cord

## **Adjustment of the water flow**

The system allows control of the flow of coolant water to each handpiece. Use the controller to complete the adjustments:

1. Lift the handpiece from the support and pull the cord forward.
2. Locate the regulator of the water flow related to the handpiece
3. Activate the water switch on the foot control.
4. Press the pedal to activate the handpiece.
5. Adjust the water flow by turning the controller clockwise to decrease the flow and counter-clockwise to increase the flow.

## **Adjustment of the air flow**

The air flow controller adjusts the air flow rate in access to the handpieces. The air flow produces the rotation of the handpieces. Use the controller to control the adjustment.

### **To adjust the flow and therefore the power of the handpieces :**

To adjust the flow and therefore the power of the handpieces:

1. Lift the handpiece from the support and pull the cord forward
2. Locate the regulator of the air flow related to the handpiece
3. Press the pedal to activate the handpiece.
4. Adjust the air flow by turning the controller clockwise to decrease the flow and counter-clockwise to increase the flow.

Through the pressure regulator, located on the side of the main unit, it is possible to check the air flow rate that is being generated. The function is to monitor the proper use of the instruments; to avoid a rapid wear of the rotors, make sure that the level of the pressure does not exceed the green line on the indicator.

## **6.2 \_ The Cuspidor**

### **6.2.1 \_ Adjustment of the Cuspidor Rinse and Filling of the Glass**

The changes of water flow to the basin and to the glass must be made within the water unit. To be able to modify the water flow, act on the valves located on the tubes that provide water to the basin or the glass. Turn the appropriate grey microball valve located on the tubes.


## 7 \_ Cleaning and Disinfection

### 7.1 \_ Cleaning and Care on the Cords of the Handpieces.

The cleaning operation of the cords can be done singularly or for all the cords simultaneously. Use the internal cleaning system powered by the specific disinfection tank. To activate the system, see the relative section in the previous pages. We suggest to repeat the operation of sanitization at least once a week; this operation will ensure the elimination of the bio-film that is deposited inside the cords and therefore a longer life of the instruments.

After each patient:

1. Remove the handpieces from the cords.
2. Fill the tank with the sanitizing solution and activate it.
3. Bring together the cords and keep them into a sink or basin.
4. Activate the cords until the tank is completely empty.
5. Repeat the operation by loading the tank with water only to rinse the tubes. Once finished the operation, turn off the tank through the appropriate switch.

Let all the air and water sources work for 20-30 seconds between a patient and another, in order to significantly reduce the biological contamination.

## 7.2 Cleaning of the Suction Cannulas

To remove the impurities that are deposited inside the suction cannulas, we suggest to aspirate two cups of cold water every two hours during the work. To guarantee the elimination of all bacteria and maintain sanitized the suction cannulas, however, we suggest to proceed to a more targeted sanitization at least once a day, at the end of all the operations on the patients.


There are specific products that are suitable for this purpose. Check on Simple & Smart's website which may be the most appropriate. In order to proceed to the sanitization procedure, perform the following steps:

1. Dilute the sanitizing solution with water into a bucket
2. Insert the suction cannulas into the bucket and hold them until the complete extraction of the solution.

## 7.3 \_ Oil Collector

Clean the oil collector (small glass placed under the main unit) once a week. For the cleaning:

1. Unscrew the glass, placed under the main unit, counter-clockwise and remove the used gauze.
2. Remove the inner sponge and rinse it using hot water.
3. Place the sponge back and screw the cover.


## 7.4 \_ Cleaning of the Suction Filter

The suction filters prevent that solids enter the suction system. A proper cleaning reduces the risk of biological contamination. It is advisable, therefore, to clean the filters daily.

To carry out the cleaning:

1. Remove each filter cap located on the outside of the water unit.
2. Pull out the filter and rinse it with warm water and a disinfectant detergent
3. Place back the filter and its cap.

Figure 19 – Suction filter

## 7.5 \_ Cleaning of the Cuspidor

The shaped spouts and the smooth surface of the basin allow a quick and easy cleaning. Remember to empty and clean the filter of the basin each time you clean the cuspidor, taking care not to empty the content inside the drainpipe.

Do not remove the spouts while cleaning the cuspidor. This will prevent that the disinfectant products, used for cleaning, cause damages to gaskets and contaminate the water for the filling of the glass.

Figure 20 – Glass basin

At the end of each day, flush the drainpipe to remove residues from the pipe. If the drain is not cleaned regularly, debris may clog the basin. To wash the drain, proceed to the rinsing of the basin several times, long enough that the water flows through the drain for at least 60 seconds. Raise the chair to facilitate the drop of the water flow.

## 7.6 \_ Cleaning of the lamp and the Unit Surfaces

In order to maintain a proper sanitization of the chair and the operating lamp, the sanitizing of all surfaces is recommended. Pay more attention to the leather seat and the glass of the operating lamp; parts which are particularly sensitive to sanitizing products in general. We suggest, on these surfaces, to simply use a cloth moistened with warm water.

For all other surfaces of the dental unit, you still need to use alcohol and ammonia free products, or minimally invasive products. Refer to the Simple & Smart's website to see which products are most suitable for such sanitization.

To clean the outer surfaces, avoid the use of aggressive products, such as Lysol, alcohol, acids, ammonia, etc.

Do not clean the glass of the lamp with product containing alcohol. Also avoid to spray liquid products directly on the buttons for the control of the movement located on the main unit and on the assistant unit.

**It is not possible to autoclave** the suction filter, the basin filter, the terminal parts of the suction cannulas and the intraoral camera.

Labels and the instructions sticked on the device have to be present and perfectly legible, if it is necessary, they will have to be completed or replaced.

**For the cleaning of the dental unit, the medical director is invited to carefully follow the instructions provided in this manual and in compliance to the laws in force.**

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Issued to

### **SIMPLE & SMART SRL**

Head office and Operative Site:  
Via Perola 16, 24021 Albino (BG)- ITALY

*Bureau Veritas Italia S.p.A. certifies that the Production Quality Assurance System of the above organization has been audited and found to be in accordance with the requirements of*

## **DIRECTIVE 93/42/EEC as amended**

*(according to Annex V)*

*In relation to the following products*

Subcategory:	Active dental devices
Generic group:	Dental units
Model:	SS_ONE, SS_KISS
Class:	Ila

Reference BV practice: ZIG. N. 60434385

Certification cycle start date: **03 July 2015**

Subject to the continued satisfactory operation of the organization's Management System, this certificate expires on: **02 July 2018**

Original certification date: **03 July 2015**

Certificate N. **IT261641**

Rev. N 2 of: 28 August 2017

**CARMELA DE FEO- Technical Committee Coordinator**

*This certificate is issued by Bureau Veritas Italia S.p.A. Via Miramare, 15 - 20126 Milan, as a notified body for the Directive 93/42/EEC, with identification number 1370*

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Declare that the product: **SS\_ONE**

**Serial Number:**  
**Production date:**

**Classification:** CLASS IIA REGULATION 5 - REGULATION 9  
**in accordance with the Annex IX DIRECTIVE 93/42 EEC. E S.M.I.**  
**conformity assessment** ANNEX V DIRECTIVE 93/42 EEC AMENDED BY  
DIRECTIVE 07/47 EC

**is accordant to the essential requirement of D.Lgs. 46/1997 about the Medical Devices Directive 93/42/CEE e s.m.i. including the directive 2007/47/CE of the D.lgs. 37/2010**

Ref CE certificate n°: IT261641  
Released by: BOREAU VERITAS ITALIA S.P.A. - Via Miramare, 15 - 20126 MILANO  
Identification number: CE 1370  
Original Certification Date: 07-03-2015

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A handwritten signature in blue ink, appearing to be 'C. Veritas', is written over the contact information.

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### **WARRANTY**

The SIMPLE & SMART equipment are guaranteed against defects in material and workmanship for a period of 12 months from the time of installation. Any warranty extensions are allowed only if explicitly written in the contract of sale or invoice. The warranty is on the dental unit, that is the dental chair excluding: handpieces , electric micromotor, Endoscopy System (camera and monitor). The warranty covers only those parts which are not subject to damages due to the use. Any costs related to shipping spare parts, even if under warranty, remain charged to the customer. The complaint of the defects on the equipment must be done within 15 days from installation.

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### **CONTROL CERTIFICATION**

SIMPLE & SMART srl certifies that the dental units SS\_ONE and SS\_KISS, in addition to the testing carried out during manufacturing, were subject to a further quality control at the logistics center, following a strict quality control process, AS APPROVED BY THE CERTIFICATE OF CONFORMITY OF THE MACHINE. This certificate can be requested at your discretion.

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**ELECTRICAL TEST**

Electrical inspections must be carried by specialized technician and are borne by the local distributor of SIMPLE & SMART products. Electrical tests must be carried out periodically, according to the current regulations and must be made with a special tested and certified tool. SIMPLE & SMART disclaims any liability if the tests are not carried out, or if are not compliant to the the legislation.

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## **ASSISTANCE**

The assistance service must be guaranteed by the dealer from whom you purchased the equipment, if the machine is not purchased directly from the manufacturer. However, to receive a more complete support service, you can also contact the service center of SIMPLE & SMART by email: [assistance@simplesmart.it](mailto:assistance@simplesmart.it) or by telephone. More details on the service can be found directly on the site [www.simplesmart.it](http://www.simplesmart.it)

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### **INSTALLATION SHEET**

In the following pages, you will find the form to fill in concerning the correct installation. The enabled technician is required to fill it after installation. The form must be signed by the technician and the customer and sent electronically to one of the simple & smart mail addresses. A copy of the form must be kept by the customer. The correct installation form is necessary in order to obtain warranty service.

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Herewith, the above mentioned installer, declares to have installed and tested the machine correctly, in accordance with the procedures and technical parameters of the device itself. The customer also declares to be aware of how to use the dental unit, of the conditions contained in the following warranty and to have read the manual.

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FIRMA/Signature/Firma  
INSTALLATORE/Installer/Instalador

FIRMA/Signature/Firma  
CLIENTE/Customer/Cliente

**Simple&Smart Srl**

Via Perola 16  
24021, Albino, Bergamo.  
ITALIA  
[www.simplesmart.it](http://www.simplesmart.it)