AIR-FLOW® CLASSIC

Before using this product, please carefully read and follow the recommendations of the instruction manuals originally delivered with your device and system. Please pay special attention to the safety precautions.

Product Description

The AIR-FLOW® CLASSIC (prophylaxis powder) is used only for supragingival treatments.

This powder is for use in dentistry with the EMS air polishing devices.

Applications

This product is recommended for:

- Removing stains/colorations and dental plaque
- Plaque removal for placement of sealants,
- Surface preparation prior to bonding/cementation of inlays onlays, crowns and veneers,
- Surface preparation prior to placing composite restorations,
- Cleaning prior to bonding orthodontic brackets,
- Effective plaque and stain removal for orthodontic patients,
- Stain removal for shade determination,
- Plaque removal prior to fluoride treatment,
- Plaque and stain removal prior to whitening procedure.

Contra-indications

Patients and users suffering from chronic bronchitis or asthma must not, under any circumstances, be treated or use an air polishing device. The jet of air and powder could cause respiratory difficulties.

Patients suffering from upper respiratory tract infections.

Patients on a low salt diet (for example hypertension) must not be treated with the powder as it contains sodium bicarbonate. For these patients it is recommended to use the salt free powders also provided by EMS (see below).

The flavors contained in the powder can lead, in specific case, to allergic reactions. For patients known to suffer from such reactions, use the AIR-FLOW® CLASSIC Neutral powder provided by EMS.

Do not direct the jet of powder toward fillings, crown or bridgework as this could damage these restorations.

Do not use this powder subgingivally. For this application use only the powder recommended by EMS for subgingival application.

Recommendations for the treatment

Patient and user have to wear eye protection.

User has to wear a mask and gloves.

Use of long cotton rolls in the mouth for maximum protection.

Use the high-speed evacuator of your dental unit.

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General information

- Do not exceed the powder filling limit as indicated on the powder chamber.
- This powder is sugar free.
- If the bottle is left opened, there is a risk of agglomeration of the powder, that could block the device.

Storage and shelf-life

- The powder should not be stored near acids or heat sources.
- Storage condition -10°C to 30°C.
- Do not use after expiration date.

Also available:

- AIR-FLOW® CLASSIC available in Lemon, Mint, Cherry and Neutral (for supragingival treatments).
- AIR-FLOW® SOFT (salt free for supragingival treatments).
- AIR-FLOW® PERIO (salt free for subgingival treatments).
- AIR-FLOW® PLUS (salt free for subgingival and supragingival treatments).

Powder composition:

- Sodium bicarbonate
- Amorphous silica
- Sweetener and aroma (Not present in AIR-FLOW® CLASSIC Neutral)

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