# OBTURYS

**FR**

#### INDICATIONS

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**EN**

#### INDICATIONS

Obturys is an epoxy-amine resin based cement for permanent root canals sealing, to use in combination with gutta-percha points. It has the following properties:

- Radiopaque
- Long term stability
- Flexibility that prevents cracking of fully cured material
- Low shrinkage
- Good flow properties

#### CONTRAINDICATIONS

Hypersensitivity to epoxy resins and amino groups. Caution: do not use in case of patient with allergic sensitivity to epoxies and amines.

#### ADMINISTRATION MODES

- Automatic mixing (Automix) : 5ml paste-paste syringe with automix tips with extra-fine Intra oral tips

- Manuel mixing (Handmix) : 5ml paste-paste syringe with specific tip for manuel mixing

#### PROUCT COMPOSITION

Paste A - base : DGEBA ; Zirconium oxide ; Bismuth oxychloride ; Ytberium oxide ; Fumed silica

Paste B - catalyst : Diamine ; Zirconium oxide ; Bismuth oxychloride ; Ytterbium oxide ; Fumed silica ; Silane quaternary ammonium salt

#### INSTRUCTIONS FOR USE

##### • Automix

Prepare the root canal, clean and dry before use of Obturys.

Remove the cap by turning it 90° counter clockwise and pulling it off. Prior to use, bleed the syringe.

Attach the mixing tip to the syringe and then turn the tip 90° clockwise. Automix syringe ensures proper homogeneous mixing of two pastes and allows direct application of the sealer into the root canal.

After the application, the mixing tip must be removed by rotating it 90° clockwise and pulling it off. Then the Automix syringe must be closed with the original cap.

##### • Handmix

Prepare the root canal, clean and dry before use of Obturys.

Remove the cap and extrude the same quantity of the two pastes on the mixing pad provided for this purpose. Then mix the two pastes using spatulas. After use, close the Handmix syringe with its original cap.

#### APPLICATION

As a standard technique, Obturys is used in combination with gutta-percha points. For filling techniques where most of the canal is obturated by endodontic point material, apply only a light coating of Obturys to the canal walls. For the Master Point-Technique, select a gutta-percha point (or alternatively a paper point or a reamer) in the size of the last instrument used during apical preparation. Wet the canal walls with Obturys through a pumping or simultaneously rotating movement in a counter clockwise direction of the point/reamer. Alternatively, apply Obturys onto the tip of a Lentulo spiral. Advance the Lentulo spiral slowly to the apex running at very low speed. Avoid the formation of bubbles in the material and overfilling of the canal. Withdraw the Lentulo very slowly, still running at low speed. Dip point into Obturys and insert them into the canal with a slow pumping motion.

#### REMOVAL OF OBTURYS

The conventional technique of removing the gutta-percha is applicable when Obturys is used with gutta-percha points.

#### PRECAUTIONS

- Under certain storage conditions, one paste of Obturys may slightly separate. It has been shown that this process does not negatively affect the performance of the mixed product.

- Inadequate cleaning may slow epoxy polymerization (setting process). Please make sure to wash the canal and dry well prior to use of Obturys.

- You may observe some brown discoloration in one paste of Obturys. This is normal behaviour of the material and has no effect on its mechanical properties.

- After mixing the discoloration will disappear and blend into the dentin shade color of the mixed material.

#### STORAGE

Obturys must be stored at temperature from 4 to 24°C (39 to 75°F). Automix and Handmix syringes must be closed after use only with original cap. Proper storage conditions guarantee 2 years shelf life of Obturys.

#### SETTING TIME

The setting time of Obturys is 48 hours at 37°C.

#### CLEANING OF INSTRUMENTS

All the instruments that were in contact with Obturys should be cleaned immediately after use with alcohol or acetone.

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