Chemin du Verger 3 – CH 1338 Ballaigues
Switzerland
www.dentsplymaillefer.com

Product : **GLYDE™ DOSES**

Ref doc. : F19 02 41.X / 05 / 2003 – updated 04/2008

CE
0459

		Page
F	PROTOCOLE D'UTILISATION	2
GB	DIRECTION FOR USE	4
D	GEBRAUCHSANWEISUNG	6
E	INSTRUCCIONES DE USO	8
I	ISTRUZIONI PER L'USO	10
NL	GEBRUIKSVOORSCHRIFT	12
DK	BRUGSANVISNING	14
FIN	KÄYTTÖOHJEET	16
S	BRUKSANVISNING	18
P	PROTOCOLO DE UTILIZAÇÃO	20
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΓΙΑ	22

GLYDE™ DOSES

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F

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DIRECTION FOR USE Root Canal Conditioner

0) COMPOSITION:

Glyde™ Root Canal Conditioner is a combination of Ethylenediametetracetic Acid (EDTA) and carbamide peroxide in a water soluble base.

The EDTA chelates calcium salts from calcified areas within the root canal to facilitate the use of files.

1) INDICATIONS FOR USE:

These products are to be used only in a clinical or hospital environment, by qualified users.

Application field:

Glyde™ is a device in the form of a gel for use in the chemical and mechanical cleansing of the root canal preparation during endodontic therapy.

Glyde™ is a combination of EDTA and carbamide peroxide specifically formulated to provide cleansing of the root canal preparation. The polyethylene oxide base used in Glyde™ acts as a lubricant to facilitate shaping of the root canal using rotary and hand instrumentation.

Glyde™ allows for cleansing action that facilitates easy removal of vital pulp tissue and necrotic pulp tissue. Glyde™ is designed to be used with endodontic irrigation with sodium hypochlorite solutions. Oxygen bubbling occurs through the release of oxygen from the carbamide peroxide. This action allows for pulp tissue, dentinal shavings and debris to be easily removed.

Glyde™ encourages lightening of the tooth if discoloration exists from non vitality. The use of sodium hypochlorite with Glyde™ promotes internal bleaching of the tooth. This process is enhanced by the release of oxygen from the carbamide peroxide.

Also, Glyde™ leaves an appropriate internal canal surface for penetration of other medicaments used in endodontic therapy by more completely cleaning the canal and opening dentinal tubules.

2) CONTRA-INDICATIONS:

Do not use Glyde™ outside of the root canal system or on soft tissues.

Immediately wash skin with soap and water after contact or flush eyes if contact occurs.

3) WARNINGS:

1. Not injecting material under pressure into the canal.
2. Carbamide Peroxide and EDTA should not be applied to tissue outside of the Root Canal such as mucosa of skin. Wash skin with soap and water after contact.
3. Glyde™ canal conditioner should extrude easily; DO NOT USE EXCESSIVE FORCE.
4. A single-unit dose of Glyde™ is designed to be used during approximately one endodontic case (tooth) with multiple canals. Use one single unit dose of Glyde™ per patient.
5. Do not take internally.

4) PRECAUTIONS:

The viscosity of Glyde™ is affected by storage conditions. If a more viscous or thicker consistency is desired, store in a refrigerator.

In the event of ingestion, contact with eyes or skin, please contact a doctor immediately.

It is recommended to wear protective clothing when using this product.

5) ADVERSE REACTIONS:

- Irritating to eyes (R36) – Irritating to skin (R38)
- Vapours may cause drowsiness and dizziness (R67)

6) STEP BY STEP INSTRUCTIONS:

1. Open the pulp chamber in the usual manner and obtain adequate access to the root canal or canals.
2. Remove the pulp and flush canal with sodium hypochlorite solution.
3. Dispense Glyde™ either directly from the unit dose applicator in and around the canal orifice,
OR
Dispense Glyde™ from the applicator to a dappen dish and carry to the canal by coating a small hand file.
4. The sodium hypochlorite will react with the Glyde™ and begin to effervesce or bubble. PLEASE TAKE NOTE OF THE BUBBLING REACTION.
5. Using a small hand file, explore the canal, mixing the Glyde™ and sodium hypochlorite.
6. Gently irrigate and evacuate with sodium hypochlorite.
7. Coat the first file in the rotary file technique with Glyde™ and use to desired depth. Glyde™ will now act as a file lubricant while instrumenting the canal with a rotary instrument.
8. Repeat step 7 for each step in the crown-down technique sequence until canal shaping and cleaning is completed. NOTE: irrigate and evacuate using sodium hypochlorite to remove Glyde™ mixture and dentinal debris from canal after every 2 files in the rotary file sequence.
9. Once shaping is completed, thoroughly flush and evacuate canal using H₂O to remove all irrigants. Dry canal and obturate.

STORAGE AND SHELF LIFE:

Glyde™ must be stored at a temperature between 2 and 12°C. Do not expose Glyde™ to temperatures exceeding 50°C.

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