**VOCO****Rebilda DC**

In accordance with DIN EN ISO 4049

## Instructions for use

## Product description:

**Rebilda DC** is a dual-curing, highly radiopaque flowable composite with excellent mechanical properties for core-build up.

**Rebilda DC** is available in three different shades to cover a large range of indications:

- blue to visualise the transition between material and tooth substance, especially in molar and pre-molar areas
- dentine shade for high aesthetic demands, e.g. under full ceramic restorations or composite restorations with high translucency
- white as a possibility to visualise preparation margins and simultaneously meet aesthetic demands

## Indications:

- Adhesive core build-up of vital and non-vital teeth
- Luting of fibre-reinforced resin posts

## Method of use:

- A) Take appropriate measures in areas close to the pulp.

## B) Bonding material

**Rebilda DC** is used in the adhesive technique in combination with a light- or dual-curing adhesive which must be compatible with dual-curing composites. The product's instructions for use should be consulted to determine the correct preparation (etching technique) and procedure.

## C1) Rebilda DC in the cartridge

Insert **Rebilda DC** cartridge into the **VOCO** dispenser (type 3) or any other extruder gun with a suitable piston. Due to technical reasons slight differences in the filling level between the cartridge chambers can occur prior to the first use. Therefore remove cartridge cap and squeeze the material out until the material is delivered evenly from both chambers. Then insert a mixing tip type 12 and lock with a 90° turn clockwise. If necessary, attach an additional intraoral tip type 3 onto the mixing tip. For application in the root canal, use a mixing tip type 19 with integrated metal cannula. The cannula can be carefully pre-bent by hand. By pressing the lever the material is extruded and automatically mixed in the correct ratio.

## C2) Rebilda DC in the QuickMix syringe

Due to technical reasons slight differences in the filling level between the syringe chambers can occur prior to the first use. Therefore remove syringe cap and squeeze the material out until the material is delivered evenly from both chambers. Afterwards place a mixing tip type 11. Match the nooks on the syringe and the mixing tip and fix the mixing tip with a 90° turn clockwise. If necessary, put an application tip type 4 on top of the mixing tip. The material is mixed automatically in the correct ratio when dispensing and can be applied directly. All QM syringes come with a detachable piston. The piston can simply be pulled off from the QM syringe and reinserted.

The plugs in the syringe chambers will stay in place preventing any cross contamination of the material in the syringe.

## D) Application of the core build-up material

Place a matrix around the prepared tooth and apply **Rebilda DC** directly from the cartridge/QuickMix syringe.

## Light-curing

The use of transparent matrices is advantageous. Place the tip of the polymerisation lamp as close as possible to the surface of the core build-up and polymerise with halogen light for 40 s.

If the distance is more than 5 mm, the curing depth may be reduced. In case of a larger core build-up, generally apply in layers of 2 mm and cure each layer for 40 s. After light-curing, the material can be finished immediately.

## Chemical curing

Without light-curing, the build-up can be finished 5 min after its application in the mouth.

**Rebilda DC** has a working time of approx. 2 min.

In case a preshaped core form is used, apply **Rebilda DC** directly into the preshaped core form. Then seat the crown on the prepared tooth, adjust its position and light-cure the margins for 40 s to fix the piece.

After use leave the mixing tip on the cartridge/QuickMix syringe or re-attach the cartridge cap.

## Luting of fibre-reinforced resin posts (in combination with a dual-curing adhesive):

1. When luting root posts, apply the dual-curing adhesive according to its instructions for use.

## 2. Using the cartridge:

When using mixing tip type 19, apply **Rebilda DC** directly into the root canal with the bendable metal cannula of the mixing tip. If necessary, the metal cannula of the mixing tip can be pre-bent by hand for this purpose.

When not using mixing tips type 19, proceed as follows:

After preparing the post according to the manufacturer's instructions, wet post with **Rebilda DC**. Do not apply **Rebilda DC** directly into the root canal nor use a lentalo! Insert the post making sure to leave a slight excess.

## Using QuickMix syringe:

It is recommended to apply **Rebilda DC** directly into the root canal using thin intraoral tips type 1. Do not use a lentalo! After preparing the root pin according to the manufacturer's instructions, it can be additionally moistened with **Rebilda DC**. Insert the post making sure to leave a slight excess.

3. To fix the post and cure the dual-curing adhesive applied to the occlusal surface, light-cure for 40 s.

4. For core build-up, follow the instructions given in D) Application of the core build-up material.

## Warnings/precautionary measures:

- For complete polymerisation, avoid any contact of the adhesive and **Rebilda DC** with phenolic substances (e.g. eugenol).
- Do not use bonding agents which are not compatible with chemically curing or dual-curing composites.
- Insulate the core build-up before taking impressions or making temporary crowns/bridges (good moistening is often sufficient) to avoid bonding to these materials.
- **Rebilda DC** contains Bis-GMA, diurethane dimethacrylate, BHT and benzoylperoxide. **Rebilda DC** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Avoid contact with the oral mucosa.

The physical values of **Rebilda DC** correspond to DIN EN ISO 4049.

## Information on storage and application:

Close **Rebilda DC** cartridges/QuickMix syringes after use. Store cool. Do not expose **Rebilda DC** to direct sunlight or the chairlight. Store at temperatures of 4 °C to 23 °C. Do not use after expiry date. Do not sterilize the mixing tips in the ultrasonic bath.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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