Circuiti e maschere per
Circuit and mask for
Schlauch mit maske für
Circuitos y mascarillas para
Circuitos e máscaras para

Master Flux

Plus

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1. INTENDED USE
3. CONNECTION
4. CLEANING

1 INTENDED USE

The Tecno-gaz circuits and masks are accessories of device Master Flux and should only be used with the device Master Flux for the supply of medical gases to the patient.

TECNO-GAZ is not responsible if the circuits and masks are used in a manner different from that stated above.

TECNO-GAZ is not liable for misuse of circuits and masks.

The masks are sold separately from the MASTER FLUX ONLY AS PARTS OF THE DEVICE.

2 | **CIRCUIT AND MASK**

Components CICUIT

<p>MF820ZMF</p>	<p>CIRCUIT COMPLETE</p>	
<p>MF835ZMF</p>	<p>MASK size 1</p>	
<p>MF836ZMF</p>	<p>MASK size 2</p>	
<p>MF837ZMF</p>	<p>MASK size 3</p>	
<p>MF838ZMF</p>	<p>MASK size 4</p>	
<p>SMFA197</p>	<p>CIRCUIT EXHAUST PASSIVE</p>	
<p>MF824ZMF</p>	<p>CIRCUIT CONNECTION OF ASPIRATOR</p>	

The circuit **MF824ZMF** is composed of:


<p>TERMINAL ASPIRATION</p>	
<p>REDUCTION FOR SUCTION SALIVA</p>	
<p>ONE-WAY VALVE</p>	

the circuit **MF820ZMF** is composed of:

<p>Fitting connection to flowmeter box ▶ Connection 22 M/F 90°</p>	
<p>Two tubes length 1.8 m (material PVC/PE) ▶ Transparent and flexible tube. ▶ Corrugated tube (rigid spiral) externally (anti-crushing) smooth inside for low compliance and noise. ▶ Flexible connectors for easy insertion.</p>	
<p>Exhaust fitting connection (material: stirolox) ▶ Reduction 22-15</p>	
<p>Reducing exhaust connection (material: TPE, PP) ▶ Elastomeric connector 9-11 mm</p>	

<p>Pipe clip support and velcro band ▶ Allows correct positioning of the tubes, the adhesive clip can be secured behind the headboard or via the Velcro strip can be positioned on the headboard without a secure fixing. ▶ It allows a positioning with fixed or rotary movement, according to the needs.</p>	
--	---

The mask is composed of:


<p>Kit Silhouette nasal and Support front Code 2ZMFA0050 size 1 Code 2ZMFA0051 size 2 Code 2ZMFA0052 size 3 Code 2ZMFA0054 size 4 (material: silicon) ▶ It reduces N2O environment emissions ▶ Create a mask placement comfortable</p>	
<p>Shell (material: PC) ▶ Transparent</p>	
<p>Small support front (material: PC) ▶ Transparent For mask size 1 and size 2</p>	
<p>Large support front (material: PC) ▶ Transparent For mask size 3 and size 4</p>	
<p>Fitting Y (PC) with valves (PC/silicone) ▶ Allows the correct flow of the gas mixture by preventing the exhaled gas can flow back into the supply circuit. ▶ The small room eliminates the risk of permanence CO2</p>	
<p>HEAD RIBBON WITH FASTENER FOR HEAD RIBBON ▶ The mask stabilizing and regulating via a velcro the length to adapt it to the patient.</p>	

CONNECTION INSTRUCTIONS

Check the integrity of the components.

Choose the most suitable mask for the patient.

Place the mask and adjust the head ribbon. The mask must be stable and comfortable for the patient. Connect the pipes, paying attention to the arrow on the Y fitting. The arrow indicates the flow direction. The non-return valves, placed in the fitting, push the gas flow in the direction indicated by the arrow.



Gas inlet pipe: Connect a hose to the special corrugated tube (Item G) to the flow control box by inserting the connector connection flow control box ".

Connect the free end to the Y connector placed on the mask, WITHOUT SIDE ARROW.

Gas inlet pipe: Connect the other corrugated pipe to the appropriate hose (Pos P) through the elastomeric connector. Connect the free end to the Y connector placed on the mask, SIDE WITH ARROW.

Fix the clip to stabilize the pipes behind the headboard of the chair.

And 'possible fixing with adhesive (fixed) or through the appropriate supplied Velcro (detachable).

EVACUATION EXHALED GASES

The operator can choose the best solution for processing the exhaled gases, using the following connections:

PASSIVE EXHAUST	
UNIT-MOUNTED MASTER FLUX	WALL-MOUNTED MASTER FLUX
<p>Connect the hose, exhaust gas circuit mask to the hose in the center at the center of the unit.</p> <p>Connect the circuit SMA197, to the hose connection at the bottom of the cabinet.</p> <p>Connect the end of the spiral tube to the outside and make sure there is a final filter.</p> <p>The discharge must take place in an open, well-ventilated.</p>	<p>Connect the exhaust gas pipe of the appropriate circuit hose (SMFA257) after the spiral tube (SMFA197).</p> <p>Connect the end of the spiral tube to the outside and make sure there is a final filter.</p> <p>The discharge must take place in an open, well-ventilated.</p>
	


ACTIVE EXHAUST

UNIT-MOUNTED MASTER FLUX

Connect the hose, exhaust gas circuit mask to the hose in the center of the cabinet.
 Connect the circuit MF824ZMF to the hose on mobile.
 Connect the terminal to the suction.

WALL-MOUNTED MASTER FLUX

MF824ZMF connect the circuit to the discharge circuit.
 Connect the terminal to the suction.

SETTING THE VALVE

Turn the ring perforated valve (SMFA236) until closing holes.
 Start dispensing OXYGEN setting ~ 10 NI / min.
 Start the surgical aspirator, set the minimum intake.
 Observe THE BAG ANESTHESIA.
 Turn the ring perforated valve gradually opening the holes until you get the ideal state of the BAG ANESTHESIA "HALF FULL HALF EMPTY".

CAUTION:

By varying FLOW, it may be necessary, adjust the correct intake flows of exhaled gases.

4 CLEANING**Use detergents that are not aggressive to plastic**

Clean the surfaces of the device with a dry cloth.

Before and after use, clean, clean or sterilize the parts in contact with the patient or with his breath.

Cleaning table for MASTER FLUX parts

			STERILISATION	DISINFECTION
CODE	DESCRIPTION	MATERIAL	STEAM	QUATERNARY SALTS
CM84001	BAG ANESTHESIA	RUBBER (NO LATEX FREE)	NO	SI
MF820ZMF	CIRCUIT	PVC	NO	NO
	MASK	PC+SILICON	T134	SI
	HEAD RIBBON	VELCRO	NO	NO
SMFA236	ONE-WAY VALVE	AL + SILICON	T134	SI
SMFA197	INTAKE	PVC + AL	NO	NO

MASK ASSEMBLY

After cleaning the device, assemble the mask according to the diagram below.

**ATTENTION TO THE DIRECTION OF FLOW
ARE THESE VALVES
BINDING THE FLOW IN THE DIRECTION SHOWN**

Circuit INLET

INLET

Circuit OUTLET

OUTLET

Exhaled Gas Circuit identification

For easy identification of discharge circuit, during are provided template-circuit, assembly, two blue rings (Cod. 1ZMFA0082).

Figures A B C D show how to use the rings. The indicator is autoclavable.

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TECNO-GAZ S.p.A.