



TePe Universal Care™

Handle Polypropylene (ISCC certified).

Filaments Polyamide (PA).

Pack *Consumer pack; 1 pc/blister pack
Retail carton; 14 pcs/carton*

The blister is made of polyethylene terephthalate, (PET), partly recycled.
The backing card is made of FSC certified cardboard.
The tray in the retail box is made of corrugated cardboard.

Retail and transport cartons; Chlorine-free cardboard respectively corrugated cardboard. All from recycled material, except for the retail cartons which are from partly recycled material and partly from virgin FSC certified material.

Recycling Toothbrush in household waste.

The blister pack should be separated and recycled; the blister in container for plastic recycling and the backing card in container for recycling of cardboard.

The cartons in container for recycling of cardboard/corrugated cardboard.

Product features

- A unique angle provides easy access to the lingual and palatal surfaces of teeth and implants.
- Facilitates cleaning in narrow and hard-to-reach areas.
- End-rounded filaments to ensure safe and gentle cleaning.
- Neck can be angled for easier reach.

Examples of areas of use The slim brush head facilitates cleaning in very narrow areas, such as lingual surfaces of natural teeth, implants, or prosthetic constructions, fixed lingual braces or retainers, along the gum line, and retroclined teeth.

The brush is also intended for cleaning of lingual surfaces if reduced mouth opening capacity occurs. Additional areas of use can occur upon recommendation from a dental professional.