



TePe GOOD Tongue Cleaner

Handle	Biobased high-density poly-ethylene, HDPE
Autoclave	Autoclaving not possible. Disinfectant 95°C possible.
Pack	1 pc/blister pack, 14 blister packs/SRP. 1 colour (green) The plastic front is made of recycled polyethylene terephthalate, (PET). The backing card is made of cardboard. The retail box is made of chlorine-free cardboard. The transportation box is made of chlorine-free corrugated cardboard.
Recycling	The tongue cleaner in the usual household waste. The blister pack should be separated and recycled: the plastic front in container for plastic waste and the backing card in container for recycling of cardboard. The retail box should be placed in container for cardboard. The transportation box should be placed in container for corrugated cardboard.
Product features	<ul style="list-style-type: none">- Helps to improve bad breath and oral health by removing bacteria from the tongue.- Slim design for easy reach of the back of the tongue.- Three cleaning blades for a "triple effect" usually makes one stroke enough to clean the tongue efficiently.- User-friendly handle.- Dishwasher safe.
User tip	Place the tongue cleaner as far back as possible, with the raised blades towards the tongue and move the tongue cleaner forward.