



TePe Compact Tuft™

Handle Polypropylene (ISCC certified).

Filaments Polyamide (PA).
The tuft contains approximately 700 filaments.

Pack *Consumer pack; 1 pc/bag*
Retail carton; 25 pcs/carton
The bag is made of Polypropylene (PP).

Consumer pack; 1 pc/blister pack
Retail carton; 14 pcs/carton
The blister is made of polyethylene terephthalate, (PET), partly recycled.
The backing card is made of FSC certified cardboard.
The tray in the retail box is made of corrugated cardboard.

Retail and transport cartons for all packaging above; Chlorine-free cardboard respectively corrugated cardboard. All from recycled material, except for the retail cartons which are from partly recycled material and partly from virgin FSC certified material.

Recycling Toothbrush in household waste.

Bag in container for plastic recycling.

The blister pack should be separated and recycled; the blister in container for plastic recycling and the backing card in container for recycling of cardboard.

The cartons in container for recycling of cardboard/corrugated cardboard.

Product features

- Optimal for precision cleaning.
- Dome-shaped brush head with dense, rounded tuft.
- Ideal for cleaning along the gumline, or on the occlusal surfaces of emerging teeth or impacted wisdom teeth.
- Neck can be angled for easier reach.

Examples of areas of use

The brush is suitable for cleaning of bars, implant-supported constructions, attachments for overdentures, fixed braces, along the gum line, crowded teeth, emerging teeth, or impacted wisdom teeth. It is also ideal for brushing children's emerging molars. Additional areas of use can occur upon recommendation from a dental professional.

Developed in collaboration with dental experts

Dr. Jiri Sedelmayer,
Faculty of Dentistry,
- Hamburg University, Germany

