



TePe Angle™

Handle	Polypropylene (ISCC certified)
Cap	Polypropylene (ISCC certified)
Filaments	Polyamide (PA)
Wire	Plastic coated stainless steel

Packaging

1-pack *Consumer pack; 1 pc/bag*
Retail carton; 80 pcs/carton
 The bag is made of Polypropylene (ISCC certified).
 The insert is made of FSC certified paper.

6-pack *Consumer pack; 6 pcs/blister pack*
Retail carton; 10 packs/carton
 The blister is made of polyethylene terephthalate (PET), partly recycled.
 The backing card is made of FSC certified cardboard.
 The tray in the retail box is made of polyethylene terephthalate (PET), partly recycled.

25-pack *Consumer pack; 25 brushes with cap/bag*
Retail carton; 10 bags/carton
 The bag is made of biobased polyethylene.

Retail and transport cartons for all packaging above; Chlorine-free cardboard respectively corrugated cardboard. All from recycled material, except for the retail cartons for the 6-packs which are from partly recycled material and partly from virgin FSC certified material.

Recycling

Brush and cap

Interdental brush in household waste.
Cap in container for plastic recycling.

Pack

Flow pack with 1 brush, and 25-pack bag: container for plastic recycling.

Insert in 1-pack: container for recycling of paper.

The 6-pack blister pack should be separated and recycled: the blister in container for plastic recycling and the backing card in container for recycling of cardboard.

The plastic tray in retail boxes with the 6-packs in container for plastic recycling.

The outer cartons in container for recycling of cardboard/corrugated cardboard.







Product features

- 6 colour coded sizes.
- Extra long grip and angled head allowing easier access to the back teeth from the inside and the outside.
- Plastic coated wire.
- User-friendly handle.
- Through sustainably sourced material, the carbon footprint is 80 % lower than for the previous Angle.

Developed in collaboration with dental experts

Prof. Rolf Attström;
Faculty of Odontology,
Center for Oral Health Sciences,
Department of Periodontology;
Malmö University, Sweden

TePe Angle™

Angle						
Colour	Pink	Orange	Red	Blue	Yellow	Green
ISO size (0-5)	0	1	2	3	4	5
Wire size (mm)	0.4	0.45	0.5	0.6	0.7	0.8