



# TePe Select™ / TePe Select™ Compact TePe Mini™

<b>Handle</b>	Polypropylene (PP).
<b>Filaments</b>	Polyamide (PA). TePe Select Medium has 1400 filaments/ TePe Select Soft has 1700 filaments/ TePe Select Extra Soft has 2500 filaments/ TePe Select Compact Medium has 1400 filaments/ TePe Select Compact Soft has 1900 filaments/ TePe Select Compact Extra Soft 2760 filaments/ TePe Mini Soft has 1320 filaments and TePe Mini Extra Soft has 1760 filaments
<b>Neck</b>	Can be angled after heating.
<b>Autoclave</b>	Program: Standard 134°C. Process duration: 55 min, sterilization duration: 7 min, sterilization pressure: 3 bar / 300 kPa.
<b>Pack</b>	1 pc/bag, 100 pcs/carton (Mini: 25 pcs/carton). Mixed colours. The bag is made of Polypropylene (PP). Select and Mini with illustrations: 25 pcs/carton. 1 pc/blister pack, 12 blister packs/carton. Mixed colours. The plastic part is made of Polyethylene terphthalate (PET). The backing card is made of cardboard. The carton is made of chlorine-free corrugated cardboard.

**Recycling** The toothbrush in the usual household waste and bag in container for plastic. The blister pack should be separated and recycled: the plastic part in container for plastic waste and the backing card in container for recycling of cardboard. The carton should be placed in container for corrugated cardboard.

**Product features**

- End-rounded filaments for gentle cleaning.
- Tapered brush-head for improved access.
- Brush-head in three sizes (adult, compact and mini).
- User-friendly handle, suits both children and adults.