

ENDO RING



A small stainless steel ring with a graduated ruler and stop attached.

Designed to measure the working length of endodontic hand instruments.

The millimetre graduations are indelibly printed on using a laser.

The ruler can be swivelled to suit Finger or Thumb and any desired direction.
It can be used for right or left hand.

The ring can be adapted to any size of Finger or Thumb.

The Endo Ring can be used with all known disinfection / decontamination products and any type of sterilization.

REF	DESCRIPTION	SIZE
182110	Endo ring Set of 2	12 x 39 mm