



Product Data Sheet

BONE FILE MILLER	
 	
Partcodes and descriptions	
BF1X	Bone file #1X Miller-Colburn #522
BF21	Bone file #21 Miller #522

Instrument category	Surgical
Commercial name	Bone file
Intended use	To shape or smooth bone
Material	Stainless Steel
Sterile	No
Single use	No
Risk class	I
Packaging	Single
GMDN	31863
Legal manufacturer	Hu-Friedy Mfg.Co.LLC.
CE Mark	Yes
Declaration of conformity	Yes

REV.	Date	Version
01	27/05/2022	1st issue