


PHOTOGRAPHIC MIRRORS


Description : Double sided photographic mirrors

Dimensions :


MZP511 Ti
MZP512 Chr
MZP516 Pur



MZP541
MZP542
MZP546



MZP521
MZP522
MZP526




MZP551
MZP552
MZP556



MZP531
MZP532
MZP536




MZP581
MZP582
MZP586



MZP561
MZP562
MZP566



MZP571
MZP572
MZP576



GLASS:

Float glass
Thickness : 3mm
weight : 0,050 Kg

COATINGS :

CHROME COATING : through evaporation under vacuum
TITANIUM COATING : through evaporation under vacuum
PURE REFLECT COATING: Niobium/tantalum through evaporation under vacuum

CLASS : I

PACKAGING : plastic sleeve unit with protective foam

Mirror glass	Dimensions
Double sided float glass - Chrome, titanium or Pure reflect through evaporation under vacuum	Thickness 3 mm - weight: 0,050 kg

Marking

Identification

Manufacturer *Prodont-Holliger*
Item code
Lot number

Marking:

CE

Cleaning

All means
except metallic or abrasive
(Ex. Scotch brite)

Sterilisation

all (autoclave, ...)
except corrosive products
follow cleaning products manufacturers
instructions

Warranty

exchange if defect, except normal wear

Storage

Dry, clean place

Shipping

48 hours on stock

Maintenance

N/A

Medical product surveillance resource

Audrey Vialle
audrey.vialle@acteongroup.com

Last revised version: 2018/06