


Blunt Needles

MIRAJECT

For a Precise Application and Rinsing



A MIRAJECT ENDOTEC LUER

Special cannula for rinsing

- Especially gentle, effective rinsing procedure
- Side outlet avoids pressure to apex
- Round needle end – low risk of perforation
- Ideal for a gentle rinsing of periodontal pockets
- Individual, sterile packing for syringes with Luer tip

B MIRAJECT ENDO LUER

Blunt, bendable cannula for endodontics

- Application of liquids and pastes into the canal in finest, controllable form
- Blunt and bendable
- Easily to insert into the root canal
- Individual, sterile packing for syringes with Luer tip

C MIRAJECT P / P SUPER D

Cannula for application of paste-like materials

- Blunt and bendable
- For cylinder ampoules (ampoule syringe), e.g. hydrocolloid, Dontisolon, etc.
- Individual, sterile packing for syringes with Luer tip

E MIRAJECT PL SUPER

Blunt application cannula

- For application of liquids, gels and for rinsing of wide root canals
- Individual, sterile packing for syringes with Luer tip

F MIRAJECT ENDO RINSE

Convenient rinsing system for root canals, etc.

- Consisting of disposable syringe and blunt cannula
- Syringe: for 2 ml of liquid with 1/10ml calibration
- Cannula: blunt, especially fine, long and bendable
- Individual, sterile packing

CE 0197

2007062.101

Need further infos? Catalog?

Phone +49 203 99269-0

Fax +49 203 299283

E-Mail export@hagerwerken.de

Web www.hagerwerken.de

P.O.B. 100654 · 47006 Duisburg, Germany



**HÄGER
WERKEN**

Dentalprodukte weltweit
Dental Products Worldwide