

1) Opening of the interdental space with OS08OP-DS 2) Contouring with OS40M-DS 3) Controlled finishing with OS25F-DS 4) Polishing with OS15POL-DS


Clinical pictures: Dr. Francesco Garino, Turin, Italy

Ref. OS	08	80	60	40	25	15
494	OS08OP-DS					
524		OS80XC-DS				
514			OS60C-DS			
514				OS40M-DS		
514					OS25F-DS	
504						OS15POL-DS
	Opening	Reduction	Contouring	Contouring	Finishing	Polishing
Diamond-coated on both sides						

Intensiv Ortho-Strips System, Double-Sided, Ref. 060B  
Assortment of 4 oscillating strips OS60C-DS, OS40M-DS, OS25F-DS, OS15POL-DS

All Intensiv Ortho-Strips Double-Sided are available in packages of 1, 3, 6 or 12 pcs

\* Intensiv Ortho-Strips System OS08OP-DS instead of Ortho-Strips OS90



1-3) Contouring of a canine mesially with OS40M-L

Clinical pictures: Dr. Francesco Garino, Turin, Italy

Intensiv Ortho-Strips System, One-Sided, Ref. 064  
Assortment of 6 oscillating strips OS40M-L, OS25F-L, OS15POL-L, OS40M-R, OS25F-R, OS15POL-R

All Intensiv Ortho-Strips One-Sided, are available in packages of 1, 3, 6 or 12 pcs

Ref. OS	08	80	60	40	25	15
494	OS08OP-R OS08OP-L					
524		OS80XC-R OS80XC-L				
514			OS60C-R OS60C-L			
514				OS40M-R OS40M-L		
514					OS25F-R OS25F-L	
504						OS15POL-R OS15POL-L
	Opening	Reduction	Contouring	Contouring	Finishing	Polishing
Diamond-coated on one side. L = Left R = Right						