

Compact stainless steel tray with colored silicone holders for the support of Intensiv Ortho-Strips during interproximal reduction (IPR) in Orthodontics

During treatments that involve the use of Intensiv Ortho-Strips, it is important to benefit from a support which is safe, hygienic, sterilizable and customizable. The identification of the file to be used has to be clear and immediate. The storage of the files after sterilization must be carried out in complete safety.

Product description

- Tray in stainless steel, consisting of a base and lid.
- The Tray has an internal folding support with six holes suitable to support sterilisable colored silicone holders.
- Dimensions: 87x51x14 mm (LxWxH).
- The Tray has slots for the disinfection liquid to flow through it.
- The Tray is designed to hold, through its silicone holders, six Intensiv Ortho-Strips, for mechanical Orthodontic stripping.

Indications

- Arrangement of Intensiv Ortho-Strips during Orthodontic stripping procedures
- Support of Intensiv Ortho-Strips during disinfection and sterilisation processes
- Storage of Intensiv Ortho-Strips

Benefits

- Intensiv Ortho-Strips available during the treatment
- Colored silicone supports for individual combination of the needed strips
- Compact and light
- Suitable for disinfection and sterilization
- Storage in complete safety of Intensiv Ortho-Strips



- 1) Introduction of the colored silicone supports in the Tray for the positioning of the Intensiv Ortho-Strips
- 2) Selection of Intensiv Ortho-Strips according to the IPR procedure









Ref. OST400
(delivered without strips)

All colored silicone holders are available in packaging of 10 pieces or in assortment of 8 pieces.

8 silicone holders are delivered free with the Tray. The desired silicone holders are to be ordered with the Tray.

Available colored silicon holders

Bur-shank silicone holders				
Ref.	056darkblue/10	056brown/10	056red/10	056white/10
µm (the grit size is related to the strips)	80	60	40	25

Bur-shank silicone holders			
Ref.	056yellow/10	056orange/10	056ass/8
µm (the grit size is related to the strips)	15	8	