


Designation	<i>TURBINE Ti-940L</i>
Reference	3122647
Visual	
Description	Standard, high-speed turbine with ceramic bearing and fibre optic lighting: powerful, silent and even lighter.
Characteristics & technical data	<p>Turbine head: standard Speed: 280,000 to 320,000 rpm Light intensity: 25,000 LUX Ultra light: 47 g Quiet operation: ≤ 55 dB Ceramic bearings: more stable and more torque Self-cleaning head system: longer service life Number of sprays: 4 Push button Pressure: 2.5 to 3 bars Type of connection: Multiflex® Kavo®. Warranty: 12 months</p>
Storage	<u>Box containing :</u> - A user manual
Legal Information	Class IIa medical device. CE0197 marking. Certifying body : TUV RHEINLAND. European representative: Shanghai International Holding Corp. GmbH. Eiffestrasse 80, 20537 HAMBURG, Germany.