

# CERAM-ETCH

9% BUFFERED HYDROFLUORIC ACID GEL

- Reparation of fractured ceramics
- Preparation of ceramic crowns, veneers and inlay surfaces prior to cementation
- Etches ceramics of all types to create a micro-porous surface which gives a strong mechanical interlock with composite resin materials
- **Stays in place**



## Indications:

### for the intrados preparation of:

vitreous ceramics crowns, veneers and inlays before cementation with resin type composites

### for all-ceramic and metal-ceramic fractured crown reparation

## Technical Data:

9% **Hydrofluoric acid**

## Reference

### Ceram Etch

2 x 1.2 ml syringes + 4 tips.....CRAM-ETCH

Find us on [www.itena-clinical.com](http://www.itena-clinical.com)

**ITENA**<sup>®</sup>  
CLINICAL PRODUCTS

Class IIa medical device marked CE 0425. Reserved for dental professionals. Read instructions for use carefully before use.  
.itena Clinical 31 avenue Georges Clemenceau, 93420 Villepinte - FRANCE - Tel: + 33 1 45 91 61 40 - IND 1 - 06/2021