



<b>Product Designation</b>	<b>NEEDLES FOR ENDODONTIC CANAL IRRIGATION</b>
<b>Reference</b>	3125099
<b>Product picture</b>	
<b>Product description</b>	<ul style="list-style-type: none"> <li>▪ Needles for endodontic canal irrigation, individually packed, for single use</li> </ul>
<b>Characteristics and composition</b>	<ul style="list-style-type: none"> <li>▪ <b>Luer-Lock</b> fixation for a safe fixation : avoid the loosening of the needle during the rinsing, compared to the clipped needles</li> <li>▪ Garter <b>stitch</b></li> <li>▪ Efficient with a <b>rebound effect</b> : the resistance moves out the remains of the canal</li> <li>▪ Gauge/diameter: 30G – 0,30 mm</li> <li>▪ Length : 23 mm</li> </ul>
<b>Packaging</b>	<p>box of 100 needles Each needle is packed individually</p>
<b>Shelf life</b>	5 years from the manufactured date
<b>Indications of use</b>	
<b>Legal information</b>	<p>Medical Device of class I sterile under European Directive 93/42/EEC including Directive 2007/47/CE CE Marking-0197-TUV RHEINLAND Manufacturer : CDM CENTER OF EXCELLENCE LTD</p>