



Nomination	ENDO PREP							
Reference(s)	<table border="1" style="width: 100%;"> <thead> <tr> <th style="text-align: center;">Products</th> <th style="text-align: center;">References</th> </tr> </thead> <tbody> <tr> <td>ENDO PREP EDTA GEL NEW CLINIX (2X3ML)</td> <td>3130544</td> </tr> <tr> <td>ENDO PREP EDTA SOLUTION NEW CLINIX (20ML)</td> <td>3130545</td> </tr> </tbody> </table>	Products	References	ENDO PREP EDTA GEL NEW CLINIX (2X3ML)	3130544	ENDO PREP EDTA SOLUTION NEW CLINIX (20ML)	3130545	
Products	References							
ENDO PREP EDTA GEL NEW CLINIX (2X3ML)	3130544							
ENDO PREP EDTA SOLUTION NEW CLINIX (20ML)	3130545							
Visual								
Description	For the widening and chemical reaming of root canals.							
Characteristics and Technical Data	<p><u>Properties:</u></p> <ul style="list-style-type: none"> Ethylene diamine tetra acetic acid (EDTA) and its derivatives have a strong chelating power on different ions. It has no fatty or water-insoluble components: a single rinse removes it completely from the canal. Endo Prep EDTA Gel contains carbamide peroxide, which forms bubbles during conventional irrigation of the canal with hypochlorite, allowing a good cleaning of the root canals <p><u>Indications:</u></p> <ul style="list-style-type: none"> Facilitates catheterization and the shaping of narrow or difficult to access canals. Solubilizes organic debris and facilitates its elimination. Lubricates root canal preparation instruments. Bactericidal. Non-toxic <p><u>Composition:</u></p> <ul style="list-style-type: none"> EDTA Na2 15%, benzalkonium chloride 9%, neutralized excipients. Contains "sodium hydroxide" for PH adjustment 							
Conditioning	Available : -in 20 ml bottle -in gel with 2 x 3ml syringes.							
Durability	2 years							
Indications for use								
Legal information	Class IIA medical device under European Directive 93/42/EEC including Directive 2007/47/EC Marking CE0318							