



<b>Designation</b>	<b>Hydrophilic gauze swabs non sterile Medibase</b>											
<b>Reference(s)</b>	<b>3131818</b>											
<b>Product picture</b>												
<b>Product description</b>	<p>Pure cotton gauze swabs in 17 threads</p> <table border="1" data-bbox="719 1055 1441 1216"> <thead> <tr> <th>Number of threads</th> <th>Number of thickness</th> <th>Folded dimensions</th> <th>Unfolded dimensions</th> </tr> </thead> <tbody> <tr> <td>17</td> <td>8</td> <td>5 x 5 cm</td> <td>14 x 20 cm</td> </tr> </tbody> </table>				Number of threads	Number of thickness	Folded dimensions	Unfolded dimensions	17	8	5 x 5 cm	14 x 20 cm
Number of threads	Number of thickness	Folded dimensions	Unfolded dimensions									
17	8	5 x 5 cm	14 x 20 cm									
<b>Characteristics and composition</b>	<ul style="list-style-type: none"> <li>- Perfectly white without pigment</li> <li>- No waste of combing or other</li> </ul> <p>Possible to sterilize with: Vapour, Ethylene Oxide, Radiation</p>											
<b>Packaging</b>	10 packs of 100 swabs each ( 1 000 swabs)											
<b>Shelf life</b>	Not applicable											
<b>Indications of use</b>												
<b>Legal information</b>	<p>Medical Device of class I under European Directive 93/42/EEC including Directive 2007/47/CE            CE Marking            Manufacturer: SYLAMED</p>											