



<b>Designation</b>	<b>Non woven swabs, non sterile Medibase</b>											
<b>Reference(s)</b>	<b>3131819</b>											
<b>Product picture</b>												
<b>Product description</b>	<ul style="list-style-type: none"> <li>• Do not realise fibers when cleaning wounds</li> <li>• Good skin tolerance, even for sensitive skin</li> <li>• Compatible with antiseptic solutions and ether.</li> <li>• Homogeneous structure allowing a sensation of softness in contact with the skin and a low risk of fraying.</li> <li>• Absence of allergenic chemical binder</li> </ul> <table border="1" data-bbox="523 1104 1506 1245"> <thead> <tr> <th>Weight/m<sup>2</sup></th> <th>Number of thickness</th> <th>Folded dimensions</th> <th>Unfolded dimensions</th> </tr> </thead> <tbody> <tr> <td>40 g</td> <td>4</td> <td>5 x 5 cm</td> <td>10 x 10 cm</td> </tr> </tbody> </table>				Weight/m <sup>2</sup>	Number of thickness	Folded dimensions	Unfolded dimensions	40 g	4	5 x 5 cm	10 x 10 cm
Weight/m <sup>2</sup>	Number of thickness	Folded dimensions	Unfolded dimensions									
40 g	4	5 x 5 cm	10 x 10 cm									
<b>Characteristics and composition</b>	Composition : 70% viscose - 30% polyester Possible to sterilize with: Vapour, Ethylene Oxide, Radiation											
<b>Packaging</b>	Box of 1 000 swabs (10 packs of 100 swabs)											
<b>Shelf life</b>	Not applicable											
<b>Indications of use</b>												
<b>Legal information</b>	Medical Device of class I under European Directive 93/42/EEC including Directive 2007/47/CE CE Marking Manufacturer: SYLAMED											