



Product	ENDOSEALER						
Reference	<table border="1"> <thead> <tr> <th>Désignations</th> <th>Références</th> </tr> </thead> <tbody> <tr> <td>ENDOSEALER SERINGUES AUTOMIX (3X5ML)</td> <td>3134118</td> </tr> <tr> <td>ENDOSEALER EMBOUTS MELANGEURS (20) + 20 IO</td> <td>3134119</td> </tr> </tbody> </table>	Désignations	Références	ENDOSEALER SERINGUES AUTOMIX (3X5ML)	3134118	ENDOSEALER EMBOUTS MELANGEURS (20) + 20 IO	3134119
Désignations	Références						
ENDOSEALER SERINGUES AUTOMIX (3X5ML)	3134118						
ENDOSEALER EMBOUTS MELANGEURS (20) + 20 IO	3134119						
Picture							
Description	Permanent obturations of root canals combined with root canal points.						
Characteristics & Technical data	<p>ENDOSEALER is a two-paste epoxy-amine resin root canal sealer in a dual auto-mixing syringe with subsequent characteristics:</p> <ul style="list-style-type: none"> • Outstanding sealing ability • High radiopacity • Not cytotoxic • Long term stability • Moderate flexibility that prevents cracking of fully cured material • Low shrinkage • Good flow properties • Complies with ISO 6876 standard <p>Composition :</p> <p>Paste 1: • Bisphenol-A-(epichlorhydrin) and epoxy resin • Zirconium Oxide • Silica • Iron Oxide pigments.</p> <p>Paste 2: • Poly(oxyethylene, oxypropylene)diamine • Silicone oil • Silica • Zirconium Oxide • Antibacterial agent.</p>						
Storage	<ul style="list-style-type: none"> • 3 syringes x 5 ml (10 g) • To be stored at temperature from 4 to 28°C. • The syringe must be closed after use with original cap only. 						
Guarantee	2 years						
Indication of use							
Legal Information	Medical device of Class I under European directive 93/42/EEC, including the directives 2007/47/CE						