




# Technical sheet

<b>Designation</b>	<b>EASYSONIC R&amp;S</b>
<b>Reference</b>	<b>3135850</b>
<b>Product description</b>	Fast small instrument cleaner and disinfectant.
<b>Picture</b>	
<b>Characteristics and technical data</b>	<ul style="list-style-type: none"><li>- Good capacity: 3 liters.</li><li>- Lightweight: 3, 8 kg.</li><li>- Time setting: up to 30 minutes.</li><li>- Operating temperature 55°C.</li><li>- Stainless steel tank easy to clean.</li><li>- Total dimensions (L x l x P) : L 256 x l 56 x 217</li><li>- Dimensions of the tank: L 23.8 x l 35 x 98</li><li>- Dimensions of the inox box: L 210 x 120 x 60</li></ul>
<b>Shelf life</b>	Guarantee: 2 years
<b>Packaging</b>	Delivered with a basket and a cover
<b>Legal information</b>	HC

Version 3.0  
Updated 2022/06/10 by Nelly S.