

870111 870112
812500A 870113

CLIP CAP PP - SAFEWEAR® PREMIUM
Medical Device Class I

FEATURES & BENEFITS

- Spunbond non woven 100% polypropylene
- Assembly by gluing and ultrasound
- Synthetic elastic WITHOUT LATEX under a hem for a better comfort
- Non-woven flexible and very resistant
- Very breathable for a better comfort
- High extensibility - can wrap all hair.

Recommendations :

- Possible use in a pharmaceutical, industrial, agro-alimentary, medical environment ...



TECHNICAL DESCRIPTION

Product name: Clip Cap PP
SafeWear® Premium

Product type: Single use, non sterile

Material : Polypropylene

Dimensions : Length : 44,5 cm (± 0,5 cm)
Width : 48 cm (± 1 cm)

Unit weight : 3,4 g (± 10%)

Originn : France



Item	Color	Packaging
870111	White <input type="checkbox"/>	1 pack = 4 dispenser boxes of 250 pieces
812500A	White <input type="checkbox"/>	Polybag of 100 pieces
870112	Blue <input type="checkbox"/>	1 pack = 4 dispenser boxes of 250 pieces
870113	Green <input type="checkbox"/>	1 pack = 4 dispenser boxes of 250 pieces

Mechanical Information (Measurements done by KOLMI Laboratory)

TEST	Standard		Average	Unit	Tolerance	
Grammage	NF EN ISO 9073-1		15	g/m ²	± 10%	
Tensile strength	NF EN ISO 9073-3	À sec	sens M:	23	N	± 10%
		À sec	sens T:	12	N	± 10%
Elongation	NF EN ISO 9073-3	À sec	sens M:	75	%	± 50%
		À sec	sens T:	75	%	± 50%

Manufacturing Microbiological Information (by MICROSEPT laboratory)

Bioburden evaluation meets the requirements of NF EN 11737:2006 Standard
Additional microbiological controls: ASR, E.coli, Staphylococcus.

870111 870112
812500A 870113

CLIP CAP PP - SAFEWEAR® PREMIUM
 Medical Device Class I

STANDARDS & CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC
 Manufactured under ISO 13485 certified quality system.
 In conformity with the relevant harmonized standard EN ISO 14971

STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

LOGISTICAL INFORMATION



Outer case specifications

Item	Size mm	Gross weight kg	QTY/pallet
870111 870112 870113	295 x 260 x 215	3,8	26 packs (=104 cases) (4 layers of 6 packs + 1 layer of 2 packs)
812500A	800 x 400 x 245	5,1	27 cases (9 layers of 3 cases)

Inner case specifications

QTE
1 pack = 4 dispenser boxes of 250 pieces
12 polybags of 100 pieces



Codes EAN

Item	Outer case	Polybag
870111	8 718 007 908 091	-
870112	8 718 007 908 114	-
870113	8 718 007 908 138	-
812500A	8 718 007 908 015	8 718 007 908 022



1 pack = 4 dispenser boxes



Medicom Healthcare BV
 Parallelweg 80a | 3931 MT Woudenberg | The Netherlands
 cs@medicom-eu.com | +31 332 866 460

