







Product	Blue polypropylene shoes cover
Reference	3139806
Picture	
Description	Blue non-woven polypropylene shoes cover.
Characteristics & Technical data	<ul style="list-style-type: none"> - Material: non-woven - Grammage : 30g/m² - Colour: blue - Weight of a shoe :4.2 grams - Good elasticity - Dimensions: 16x40 cm - Dust-proof - Anti-fouling - Suitable for production workshop and medical place
Packaging	<ul style="list-style-type: none"> - 100 shoes cover per sachet - 10 sachets per carton - Weight: 5.5 kg - Dimensions: 40x30x32 cm
Storage	The product should be stored indoors (with no corrosive gas) with good ventilation and relative humidity no more than 80%.
Indications of use	<div style="display: flex; justify-content: space-around; align-items: center;">      </div> <ol style="list-style-type: none"> 1. Store the product in a dry and ventilated place to avoid moisture. 2. This product is to be used only once and then destroyed. 3. Avoid contact with high temperature / strong light sources. 4. Not to be used beyond the expiry date.
Shelf life of the product	3 years
Legal Information	Class I medical device under European Directive 93/42/EEC including Directive 2007/47/EC Manufacturer: CDM CENTER OF EXCELLENCE LTD