


MIXPAC™ T-Mixer™ Colibri™ plus breakable



Description

The MIXPAC™ T-Mixer Colibri plus breakable simplifies selection of the mixing tip since it can be used for various materials and clinical applications, e.g. impressions of the root canals, preparations and cavities. The integrated stainless steel cannula, which is bendable and rotatable up to 360°, allows for precise clinical application with a low risk of trauma to soft tissue due to its smooth, rounded design. Its integrated breaking zone located in the anterior part of the mixing tip, can be quickly and easily removed to facilitate application of large amounts of material.

Item no.	Item	Description	Volume	Ratio	Color
165 788	MBT 4.2-12-N11B	T-Mixe Colibri plus breakable	18-75 mL	1:1 / 2:1	Mint, natural elements

Your added value

Offer your dental customers a complete system with the latest technology	<ul style="list-style-type: none"> • Fulfills user requirement toward simplification • All-in-one solution • Less inventory required • Be a pioneer of a new innovative mixing solution • Increased product brand awareness due to differentiation with a new mixing tip
A partner you can trust	<ul style="list-style-type: none"> • Customized bag packaging on request • Support in marketing and communication to end-users

What dentists say!



Dr. Eric Van Dooren
Antwerp Belgium

“Being able to use almost the same mixing tip for the majority of clinical situations is an absolute must to minimize inventory and storage requirements in my dental practice. It’s easier to use an all-in-one mixing system without the need to attach an intra oral tip in order to receive the key advantages of accurate application, while the anterior part of the mixing tip can be quickly and easily removed which allows for very fast dispensing/application.”




Prof. Dr. Carlos Eduardo Francci
Sao Paulo Brazil

“The T-Mixer Colibri plus breakable was designed to improve upon the existing innovative mixing tip technology for easier and faster delivery of the material with the integrated advantages of the T-Mixer technology.”



T-Mixer™ Colibri™ plus breakable

Product advantages and features



360° rotatable

Smooth rounded design

Bendable integrated cannula

Allows flexible and highly precise dispensing, even in difficult to access area.

Integrated breaking zone

For quick and easy application

T-Mixer technology

MIXPAC™ Stamp

Genuine MIXPAC™ mixing tips have the MIXPAC™ brand stamped.

T-Mixer Colibri plus breakable

One color for all-in-one

COLIBRI™ plus breakable



Product advantages

Integrated smoothly rounded stainless steel cannula	<ul style="list-style-type: none"> • Bendable and 360° rotatable - Easy access to difficult areas • Constant flow of material, even when the cannula is bent • No additional intra oral tip is required • Low risk of soft tissue trauma due to smooth, rounded cannula
Predetermined breaking zone	<ul style="list-style-type: none"> • Anterior part of mixing tip is removable fast and easy • Application of large amounts of material
T-Mixer technology	<ul style="list-style-type: none"> • Less material waste up to 40% • Homogenous mixing results

Sulzer Mixpac Ltd | Ruetistrasse 7 | 9469 Haag, Switzerland
 P +41 81 772 20 00 | F +41 81 772 20 01 | mixpac@sulzer.com | www.sulzer.com/dental

Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. © Sulzer 2020

