

Polishette

mini shaped polishers

- Mini shaped polishers with tapered shank necks
- Superior precision
- More space to work, better access, increased visibility
- Three grits cover all tasks from contouring to high shine finish
- Single-use for convenience



SET REF PS710X.024



REF	reorder	PS715C	PS716C	PS715M	PS716M	PS715F	PS716F							
ISO	1 - 650 244 2 - 630 244 3 - 830 244	030 543 060	243 543 030	030 523 060	243 523 030	030 503 060	243 503 030							
	∅ × L, mm	6,5 x 5	3 x 7	6,5 x 5	3 x 7	6,5 x 5	3 x 7							
μ	grit	color	rec.	max.	μ	grit	color	rec.	max.	μ	grit	color	rec.	max.
	coarse	green	5 000 RPM	5 000 RPM	medium	red	5 000 RPM	5 000 RPM		extra-fine	purple	5 000 RPM	5 000 RPM	
	Refill Box of:	N Pressure:		Single-Use		Shank:		Medical device:						
	15 (diamond), 30 (non-diamond)	0.3-1N		disposable		RA, medical plastic		Class IIa						

