



KENDA DELUXE

2-STEP DIAMOND POLISHING SYSTEM

- Deluxe double diamond combination for the most esthetic and high shine polish on all composite restorations
- Top choice for esthetic dentistry
- No additional paste for high shine polishing needed



SET REF 0400.006A



REF	reorder	0406	0408	0409	0411	0456	0458	0459	0461
ISO	1: B58 204 ... 2: B38 204 ...	243 513 030	243 513 050	303 513 100	030 513 060	243 503 030	243 503 050	303 503 100	030 503 060
	∅ × L mm	3 x 7	5 x 10	10 x 1	6 x 7,5	3 x 7	5 x 10	10 x 1	6 x 7,5
	μ grit								
	color								
	rec.								
	max.								
	μ grit								
	color								
	rec.								
	max.								
	Refill Box of:	N Pressure:		Autoclavable:		Shank:		Medical device:	
	6, 12, 25, 100	0,3 - 2N		135°C, 275°F		RA, stainless steel		Class IIa	

