




KENDA HYBRID

2-STEP POLISHING SYSTEM

- Optimized for hybrid composites and ceromer
- Firm and robust polisher facilitates exact adjustments
- Non-diamond prepolisher and diamond high shine polisher are the most economic combination to yield highly esthetic results
- Long lasting
- Gives ultra-smooth surfaces with high lustre



SET REF 790.012A



REF reorder	0005	0006	0008	0009	0010	0105	0106	0108	0109	0110
ISO 1: 632 204 ... 2: 822 204 ...	030 524 060	243 524 030	243 524 050	303 524 100	030 524 090	030 494 060	243 494 030	243 494 050	303 494 100	030 494 090
$\varnothing \times L$ mm	6,5 x 10	3 x 7	5 x 10	10 x 1	9 x 8,5	6,5 x 10	3 x 7	5 x 10	10 x 1	9 x 8,5
μ grit	color		\circlearrowleft rec.	\circlearrowright max.	μ grit	color		\circlearrowleft rec.	\circlearrowright max.	
medium	whitish grey		7500 - 10 000 RPM	10 000 RPM	ultra fine	pale pink		5000 - 7500 RPM	10 000 RPM	
Refill Box of:	N Pressure:		Autoclavable:		Shank:		Medical device:			
6, 12, 25, 100	0,3 - 2N		135°C, 275°F		RA, stainless steel		Class IIa			

