

KENDA PLANUS

2-STEP SMOOTHING AND POLISHING SYSTEM

- Polisher for intra-oral use on acrylics
- Autoclavable
- Ideal for chairside adjustments on mouthguards, retainers, acrylic dentures and even snoring mouth pieces



SET REF 0600.012



REF reorder	0601	0608	0611	0618			
ISO 1: 638 204 ... 2: 638 204 ...	243 533 055	243 533 050	243 513 055	243 513 050			
$\varnothing \times L$ mm	5,4 x 13	5 x 10	5,4 x 13	5 x 10			
μ grit	color	rec.	max.	μ grit	color	rec.	max.
coarse	lime-green	7500 - 10 000 RPM	10 000 RPM	fine	grey	5000 - 7500 RPM	10 000 RPM
Refill Box of:	N Pressure:	Autoclavable:	Shank:	Medical device:			
6, 12, 25, 100	0,3 - 2N	135°C, 275°F	RA, stainless steel	Class IIa			

