

# KENDA PROPHY

1-STEP POLISHING SYSTEM

- For Tooth Cleaning on Natural Enamel
- Gently removes extrinsic stains without polishing paste
- Fine abrasive grit is already interspersed in the head
- Is firmer than common prophycups which allows easy control for specific removal
- Ideal for small esthetic corrections



SET REF 5100.012

20



<b>REF</b> reorder	5105	5106	5108	5110
<b>ISO</b> 1: 658 204...	030 523 060	243 523 030	243 523 050	030 523 090
$\varnothing \times L$ mm	6,5 x 10	3 x 7	5 x 10	9 x 8,5
$\mu$ grit	medium		color	blue
rec.	7500 - 10 000 RPM		max.	10 000 RPM
Refill Box of:	6, 12, 25, 100	N Pressure:	0,3 - 2N	Autoclavable:
			135°C, 275°F	Shank:
			RA, stainless steel	Medical device:
				Class IIa

