

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



TECHNICAL PRODUCT SHEET

PRODUCT	Box for disinfecting dental instruments
ITEM CODE	G11
MATERIAL	ABS
THERMAL CHARACTERISTICS	Not autoclavable, melting point 90°C
PHYSICAL AND MECHANICAL CHARACTERISTICS	Can contain about 1,2 litres of disinfectant and can disinfect up to 1,5 kgs of dental instruments. The plastic, holed basket that can be risen guarantees that the dental instruments dry off correctly. Measurements: cm 26x11 h8 Size of basket: cm 23x9
ACCESSORIES ON REQUEST G11.1	An accessory that contains burs and small dental instruments Size: cm 4,1x3,8
CHEMICAL CHARACTERISTICS	Resistant to all chemical sterilizers
SAFETY DATA SHEET	N/N
PRECAUTIONS	Be careful when using strong solutions containing alcohol
ANALYSIS OF RISKS IN REF. D.L. 46 24/2/97	N/A
CLASSIFICATION CEE 93/42	It is not a medical device

N/A = NOT APPLICABLE

N/N = NOT NECESSARY

Larident Srl

Via Lamaneigra, 12 r (Ponte Rimassa) - 16030 Tribogna Genova

Tel. +39 0185 939 481 - Fax +0185 938 595

info@larident.it - www.larident.it



SCHEDA TECNICA DEL PRODOTTO

PRODOTTO	Box disinfezione strumenti
CODICE	G11
MATERIALE	ABS
PROPRIETÀ TERMICHE	Non autoclavabile punto di fusione 90°C
PROPRIETÀ MECCANICHE, FISICHE O CARATTERISTICHE	Contiene circa 1,2 litri di disinfettante e può disinfettare fino a 1,5 kg di strumenti Il cestello sollevandosi garantisce l'asciugatura degli strumenti Misure: cm 26x11 h8 Misure utili del cestello: cm 23x9
ACCESSORI SU RICHIESTA G11.1	Accessorio per contenere frese e piccoli strumenti Misure utili: cm 4,1x3,8
PROPRIETÀ CHIMICHE	Resistente a tutti gli sterilizzanti chimici
SCHEDA DI SICUREZZA	N/N
PRECAUZIONI	Attenzione alle soluzioni fortemente alcoliche
ANALISI DEI RISCHI IN RIF. D. L.46 24/2/97	N/A
CLASSIFICAZIONE CEE 93/42	Non è un dispositivo medico

N/A = NON APPLICABILE

N/N = NON NECESSARIA