

Technische Daten

SOPIRA CARPULE Needle

SOPIRA CARPULE Needles are dental needles for dental anaesthesia for single use.

Needles may be used with all commercial syringes and syringe systems in common use in dental anaesthesia.

Sopira Carpule Needles are cold rolled stainless steel needles, coated inside with Lubricant.

100 sterile individually packed disposable dental needles for single use.

- Silicone coated for easier penetration and sterile.
- Plastic hub with thread and needle bevel orientation indicator for optimum injection direction and reduces insertion pain.
- Thin-wall SOPIRA Carpule Free Flow cannula for low injection pressure.

100 sterile, individually packed disposable dental needles for single use are packed with the instruction for use in a box. The box is covered with foil.

Only use needles if the seal on the packaging is intact. Do not bend the needles. Each needle may only be used for one single patient. Avoid contact with the used needle (risk of infection, do not manually recap!) Dispose of needles in special safe containers (puncture-resistant/lockable) which are available commercially. Needles should be collected and protected against human contact (i.e. in a suitable protective container) and the disposed of in accordance with valid regulations.

Heraeus Kulzer GmbH

i.A. 
Annegrete Wegner
Manager Regulatory Affairs/Dentistry

i.A. 
Heike Jansen
Regulatory Affairs/Dentistry