



<b>Designation</b>	<b>Sterile needles Medibase (100)</b>					
<b>Reference(s)</b>		<b>Description</b>	<b>G</b>	<b>Ø</b>	<b>Longueur</b>	<b>Reference</b>
	<b><u>Périapicales</u></b>	Extra-fine	30 G	0,30	16 mm	3114384
		Extra-fine	30 G	0,30	21 mm	3114385
		Fine	27 G	0,40	8 mm	3114386
			25G	0,50	36 mm	3114387
			27 G	0,40	21 mm	3114388
	<b><u>Intraseptales</u></b>	Fine	27 G	0,40	16 mm	3114389
	<b><u>Intraligamentaire</u></b>	Extra-fine	30 G	0,30	12 mm	3114390
<b><u>Tronculaires</u></b>	Extra-fine	30 G	0,30	25 mm	3114391	
<b>Visual</b>						
<b>Description</b>	Sterile needles for single use					
<b>Features and Technical Data</b>	<ul style="list-style-type: none"> <li>- Excellent quality, with silicone coating</li> <li>- Triple bevel for a gentle injection</li> <li>- Sterile with sterilization control ring</li> <li>- With bevel indicator</li> <li>- Sterilized with ethylene oxide</li>   <li>- Cannula: stainless steel</li> <li>- Plastic part: polypropylene (PP)</li> <li>- Lubricant: silicone</li> </ul>					
<b>Packaging</b>	Box of 100 needles					
<b>Lifetime</b>	5 years from the date of manufacture					
<b>Indications of use</b>						
<b>legal informations</b>	<p>Class IIA medical device under the European Directive 93/42 / EEC including the 2007/47 / EC directives            CE Marking 0426. ITALCERT Notify body            Manufacturer: TECNOFAR SPA</p>					

